

Xmatrix® Infinity

Automated Staining System for IHC, ISH, SS and FISH

OPERATOR'S MANUAL

Cat. No.: AS4000RX



FOR RESEARCH USE ONLY



Doc. No: 953-0035.0 Rev. L. ECO: 734-SEZ

Copyright ©2010 BioGenex Laboratories, Inc. All Rights reserved.

Customer Support: 1-800-421-4149 (USA and Canada Only)

For support outside USA and Canada, please contact your local distributor

E-mail: support@biogenex.com

48810 Kato Road,
Suite 100E & 200E,
Fremont, CA 94538
Tel : +1 (800) 421-4149
Fax : +1 (510) 824-1490
E-mail:
support@biogenex.com

 Plot No: 4, Survey No: 656/A,
Aerospace and precision
engineering industries SEZ.
Adibatla (V), Ibrahimpatnam
(M), Ranga Reddy District,
Telangana, India.
Pin: 501510

EC	REP
----	-----

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

IMPORTANT:

The Xmatrix® Automated staining system must be unpacked, inspected, and installed by an authorized BioGenex representative.

Read and understand all the safety and operation instructions contained in this manual before you attempt to use this system.

Section 2 (Safety Precautions) of this manual expounds more specifically to the safety instructions that must be followed to ensure safe operation and to maintain the Xmatrix® Automated staining system in safe condition. BioGenex Laboratories, Inc. assumes no liability for the user's failure to comply with those safety precautions.

The material in this manual is subject to change without notice. There may be revisions in new editions of this manual, please contact BioGenex Customer Service to inquire about the current edition of this manual or to purchase a copy.

Xmatrix is a registered trademark of BioGenex Laboratories, Inc.


Viton is a registered trademark of E.I. DuPont de Nemours and Company.

Tygon is a registered trademark of Saint Gobain Performance Plastics

Microsoft Windows XP is a registered trademark of Microsoft Corporation.

BioGenex Laboratories, Inc.
48810 Kato Road,
Suite 100E&200E,
Fremont, CA 94538
Tel : +1 (800) 421-4149
Fax : +1 (510) 824-1490
E-mail: support@biogenex.com
www.biogenex.com

Overview of the Xmatrx® Infinity Automated Staining System Software Version 2.5R01b



BioGenex Xmatrx® Automated Staining System is an array of reliable and efficient equipment that robotically executes any slide-based staining protocol with excellent consistency and accuracy. It precisely controls and regulates the temperature of 40 slides individually in a single run.

In this manual, the entire combination of the Xmatrx® workstation, the computer system, and the bulk fluids are referred to as “BioGenex Xmatrx® Automated Staining System”, or simply “the system”.

The mechanical mainframe that directly performs automated staining is referred to as “Xmatrx® workstation”, or simply “the workstation”, and occasionally “the instrument”.

Checklist of Xmatrx® Infinity Automated Slide Staining System with Accessories

S.No.	Description of Item	Part Number	Qty
1.	Slide Carrier	6520-30704	4
2.	Reagent 49 Vial Tray	6520-60104	1
3.	Pipette Tip Waste Box	6520-17701C	1
4.	Reagent Vial 17 ml -	6520-37836	6
5.	Cover-slip boxes/stacks (18x18mm)	6520-37981	3
6.	Cover-slip boxes/stack (25x25mm)	6520-37980	3
7.	Cover-slip boxes/stack (25x40mm)	6520-37979	3
8.	Cover-slip Disposal Box	6520-17705B	1
9.	Computer/CPU with Software Loaded	4270-00112	1
10.	Flat Monitor	4270-00113	1
11.	Printer For (IND Domestic Only)	4270-00122	1
12.	i 500 BLS-500	BLS500	1
13.	Printer Table/Stand	6520-37807	1
14.	Leveler	6520-04116	1
15.	IHC Xflex Consumable kit for Xmatrx Elite and Infinity	XT148-YCDE	1
16.	Assembly Calibration Disk	6520-30753	1
17.	Software Distribution Package for Xmatrx Infinity	4812-00115	1
18.	Xmatrx Series Instrument Installation Check List	951-6122.1	1
19.	UPS (For U.S Domestic Only)	4290-04029	1
20.	Assy Carboy 4L, IHC Buffer	6520-30715	1
21.	Assy Carboy 4L, SS Buffer	6520-30716	1
22.	Assy Carboy 4L, ISH Buffer	6520-30717	1
23.	Assy Carboy 4L, DEPC Buffer	6520-30718	1
24.	Assy Carboy 4L, DI Water	6520-30719	1
25.	Assy Carboy 4L, Alcohol	6520-30720	1
26.	Assy Carboy 4L, Dewax	6520-30721	1
27.	Acrylic Sheet	6520-20472	1
28.	Assy., Heater wiring	4501-27401	2
29.	2ml vial Fixture bottle	4911 02826	1
30.	2ml vial Graduated	4911 00024	1
31.	Fuse, 10AMP, Type MDA	4280-02245	2
32.	Operator's Manual package for Xmatrx Infinity System (Auto executable CD)	4812-00119	1
33.	Power Cord for China	1650-02985	01

34.	Power Cord, G, 10A/250V, Italy	1650-02965	01
35.	Power Cord, K, 10A/250V, UK	1650-02966	01
36.	Power Cord, C, 10A/250V, Europe	1650-02967	01
37.	Power Cord, E, 10A/250V, Australia	1650-02968	01
38.	Power Cord, A, 10A/250V, USA	1650-02970	01
39.	Power Cord, H, 10A/250V, Israrel	1650-02971	01
40.	Power Cord, B, 10A/250V, Japan	1650-02972	01
41.	Power Cord, I, 10A/250V, India/South Africa	1650-02973	01
42.	Power Cord, F, 10A/250V, Denmark	1650-02974	01
43.	Power Cord, D, 10A/250V, Switzerland	1650-02967	01
44.	Acrylic Sheet	6520-20472	01
45.	Screws 8-32, ¼" FHP	4510-08324-13	07
46.	Assay, Heater Wiring	4501-27401	02

NOTE: One line item from the following power cord list (Item# 33 to 43 , depending on the country of destination

Xmatrx Infinity								
Sl.No	Item	Catalog No.	Qty	No of Tests	IHC	ISH	SS	*FISH
1	IHC Xflex Ultra Consumable Kit contains below items)	XT148-YCDE	1 kit	200	✓	NA	NA	NA
	1a BARRIER SLIDE 25 X 40MM	XT151-SL	3	200			✓	✓
	1b COVERSLEIPS 25 X 40MM	XT118-YRK	2				✓	✓
	1c PIPETTE TIPS	XT146-01X	1				✓	✓
	1d PIPETTE TIPS	XT153-02X	2				✓	✓
2	ISH Consumable kit- Xmatrx (Contains Below items)	XT144-YAD	1 kit	100	NA	✓	NA	NA
	1a BARRIER SLIDES 25 X 25MM	XT108-SL	2	100			NA	✓
	1b COVERSLEIPS 25 X 25MM	XT122-YQK	2				NA	✓
	1c PIPETTE TIPS	XT146-01X	1				NA	✓
	1d PIPETTE TIPS	XT153-02X	2				NA	✓
3	Xviz Detection kit	QD550-YCDE	1 kit	200	✓	NA	NA	NA
4	XISH One step Polymer HRP ISH Detection kit	DF400-YADE	1 kit	100	NA	✓	NA	NA
5	Super Sensitive Wash buffer (500ml. Conc.)	HK583-5KE	1 kit	100	✓	✓	✓	✓
6	EZ Dewax/ X-Dewax	HK584-5KE / HX016-XAK	1 kit	40	✓	✓	✓	✓
7	Antibodies*	refer catalog	1 vial	50/200	✓	NA	NA	NA
8	ISH Probes*	refer catalog	1 vial	25	--	✓	NA	NA
9	Special Staining Kit*	Refer catalog for required stain	1	50	NA	NA	✓	NA
10	eFISH Histo kit / eFISH Cyto kit	DF500-20XE / DF510-20XE	1	20	NA	NA	NA	✓
11	Barrier Slides 18x18 mm	XT114-SL	1	70	NA	NA	NA	✓
12	Cover Slips 18x18mm	XT121-YBX	1	20 approx	NA	NA	NA	✓
13	DAPI	HK606-10K	1	10 tests	NA	NA	NA	✓
14	eFISH probes	Refer Catalog	1	10/20 tests	NA	NA	NA	✓
15	Optimiser Reagent vials 20ml unlabelled , White/Brown	XT026-V24 / XT101-24X	1	24 Vials	NA	NA	NA	✓
16	XMOUNT™ for Xmatrx® Infinity	HX035-10x	1	200	✓	✓	✓	NA
17	Vial Insert - 2ml 24 per pack	XT149-v24	1	NA	✓	✓	✓	✓
18	Brown Vial without neck 24 per pack	XT126-24V	1	NA	✓	✓	✓	NA
19	Unlabeled Xmatrx infinity 17ml white vial	XT154-v24	1	NA	✓	✓	✓	NA

20	Unlabeled Xmatrix infinity 17ml brown vials	XT155-v24	1	NA	✓	✓	✓	NA
21	Reagent Tray for FISH (30 Position)	6520-37404B	1	NA	NA	NA	NA	✓

NOTE: ALL THE ABOVE MENTIONED (CONSUMABLES, REAGENTS AND DETECTION KITS) ARE NOT SUPPLIED ALONG WITH THE INSTRUMENT, NEED TO ORDER AS PER THE REQUIREMENT OF THE APPLICATION.

NOTE:

- There are no user serviceable parts in the Xmatrix® system. When any service is needed, please consult BioGenex Laboratories, Inc. at 1-800-421-4149 (for USA and Canada only), or 510-824-1400 for all repairs.
- To order consumables or reagents, please contact BioGenex Laboratories, Inc. For detailed contact information, see page IV of this manual.

TABLE OF CONTENTS	
Chapter	Page No.
SECTION 1	11
INTRODUCTION TO Xmatrix®	11
I. INTENDED USE:	11
II. PRINCIPLES OF THE WORKSTATION	11
III. TRANSPORTATION AND STORAGE OF THE WORKSTATION	11
IV. SPECIFICATIONS	12
SECTION 2	13
SAFETY PRECAUTIONS	13
I. PRECAUTIONS AND WARNING	13
II. DISCLAIMER REGARDING STAINING RESULTS	15
III. PRECAUTIONS REGARDING BIOHAZARDOUS MATERIALS AND WASTE	15
SECTION 3	16
ON-SITE INSTALLATION	16
I. INSTALLATION REQUIREMENTS	16
II. INSTALLATION OF BUFFER CARBOYS TO THE LIQUID PANEL.	17
III. WASTE CARBOY CONNECTION	25
IV. GLOSSARY	27
SECTION 4	29
STARTING THE SYSTEM	29
I. LOADING CONSUMABLES	29
II. STARTING THE SYSTEM FLOWCHART	32
SECTION 5	35
MODULE A - IHC, ISH & SS Protocols	35
CONFIGURATION SETTINGS	35
I. MAINTAINING USER PROFILES	35
II. SPECIMEN NAME	37
MODULE B - FISH Protocols	39
CONFIGURATION SETTINGS	39
I. MAINTAINING USER PROFILES	39
II. SPECIMEN NAME	41
SECTION 6	43
MODULE A	43
BUILDING AND INSERTING PROTOCOLS	43
I. PROTOCOL TYPES	43
II. STAINING STEPS	43
III. BUILDING PROTOCOLS	44

IV. INSERTING A PROTOCOL	45
V. EDITING PROTOCOLS	46
MODULE B	50
BUILDING AND INSERTING PROTOCOLS	52
SECTION 7	53
MODULE A	53
USING PROTOCOL/SLIDE EDITOR WINDOWS	53
I. EDITOR WINDOWS	53
II. OPENING AND CLOSING THE <NEW PROTOCOL EDITOR> WINDOW	54
III. WORKING IN THE <NEW PROTOCOL EDITOR> WINDOW	55
IV. BUTTON OPTIONS IN THE EDITOR WINDOWS	56
V. SLIDE INFORMATION AREA	57
VI. PROTOCOL INFORMATION AREA	58
MODULE B	66
I. EDITOR WINDOWS	66
II. OPENING AND CLOSING THE <NEW PROTOCOL EDITOR> WINDOW	67
III. WORKING IN THE <NEW PROTOCOL EDITOR> WINDOW	68
IV. BUTTON OPTIONS IN THE EDITOR WINDOWS	69
V. SLIDE INFORMATION AREA	70
VI. PROTOCOL INFORMATION AREA	71
SECTION 8	72
EXECUTING A RUN	72
I. UNDERSTANDING A RUN	72
II. EXECUTING A STAINING RUN:	73
III. MONITORING A RUN	76
IV. PAUSING A RUN	77
V. STOPPING A RUN	78
VI. EMERGENCY STOP	80
VII. THE TYPICAL ERROR MESSAGES	82
SECTION 9	83
GENERATING REPORTS	83
I. LAUNCHING THE REPORTS MODULE	83
II. DEFAULT REPORT TYPES	85
III. PREVIEWING AND PRINTING A REPORT	86
IV. EXITING THE REPORTS MODULE	86
SECTION 10	87
TROUBLESHOOTING	87
I. SLIDE HEATING ERROR	87

II. ERRORS AT PIPETTE TRAY	88
III.A PIPETTE DROP ERROR (Before Reagent Aspiration)	90
III.B PIPETTE DROP ERROR (After Reagent Aspiration)	91
IV. PIPETTE DISPOSING ERROR	92
V. EMPTY COVERSIP BOX ERROR	92
VI. COVERSIP PLACEMENT / PICKUP ERROR	94
VIII. BASIC TROUBLESHOOTING	96
SECTION 11	97
ROUTINE MAINTENANCE	97
I. DAILY CLEANING	97
II. MONTHLY MAINTENANCE SCHEDULE	101
III. HALF YEARLY MAINTENANCE ARRANGEMENT	102
APPENDIX 1: SERVICE PLAN	103
APPENDIX 2: PREVENTIVE MAINTENANCE CHECKLIST FOR Xmatrix® SYSTEM	106
APPENDIX 3: SOFTWARE LICENSE	108
APPENDIX 4: IMPORTANT - READ CAREFULLY BEFORE USING THE XMATRIX® AUTOMATED STAINING SYSTEM WHICH CONTAINS MICROSOFT® SOFTWARE	110

SECTION 1

INTRODUCTION TO Xmatrix®

I. Intended Use:

Xmatrix® Automated Slide Staining System is for Research Use Only. This system performs Immunohistochemistry (IHC), *In Situ* Hybridization (ISH), Special Stains (SS) tissue staining, FISH and many more techniques.

The Xmatrix® system must be used and handled by skilled and trained personnel. All personnel who will be using the system must be trained. Failure to comply could result in damage to the instrument and injury to the user.

A qualified pathologist must interpret the staining results. This system does not provide interpretation of any staining results. Such interpretation is solely the responsibility of the user.

II. Principles of the Workstation

The workstation has an XYZ robotic arm that picks up standard disposable pipette tips, aspirates and dispenses reagents on specimens on the slide(s) placed in the workstation, and discards the pipette tips, thereby eliminating the possibility of cross-contamination.

The integrated multifunctional Z-Head rinses slide after each incubation and blows excess buffer from the slides at the end of each wash step. The system regulates the incubation time and temperature of the slide(s) and controls evaporation by covering the specimen with coverslips at any step.

The computer is for running the Xmatrix® System, and should not be used for any other purpose.

III. Transportation and Storage of the Workstation

New and disinfected workstation(s) may be transported on any commercial vehicle suitable for the purpose, with adequate protective packaging provided by the original shipping crate. If the workstation is to be stored for an extended period of time, the original shipping crate must be used. The storage environment should be indoors, dry and protected from freezing temperatures, insects and rodents.

IV. Specifications

Item	Specification	Item	Specification
Total Slide Capacity	4 removable slide carriers x 10 slides per carrier = 40 slides	Evaporation Control	Oil sealing
Reagent Capacity	1 Reagent Tray with up to 49 vials of reagents	Bulk Fluids	7 User Defined Buffers (including DeWax and Alcohol)
Reagent Dispense Volumes	10-850 µl	Staining Operation	Multiple protocols can run simultaneously Multiple parameters within a single run
Reagent Dispenser	Standard 1 ml disposable tips - for macro reagents Standard 200 µl disposable tips – for micro reagents	Weight	400 lb/182 Kg
Pipette Tip Trays	Four trays: Two 96-tip disposable trays for 1 ml tips Two 96-tip disposable trays for 200 µl tips	Electrical Specifications	115V 60Hz 900W 230V 50Hz 900W
Slide Temperature	25°C to 105°C Each slide individually temperature controlled	Ambient Operating Temperature	15-30°C (59-86°F) Humidity 15%-55% RH
Integrated Cover-slipper	Places and removes individual coverslips Two coverslip boxes each for: <ul style="list-style-type: none"> ▪ 25 x 40 mm ▪ 25 x 25 mm ▪ 18 x 18 mm 	Dimensions	Depth: 29"/73.7 cm Width: 46"/116.8 cm Height: 59"/149.9 cm
		Software	Proprietary BioGenex software running on Microsoft® Windows(XP/7/8/9/10/OT)

SECTION 2

SAFETY PRECAUTIONS

The following general precautions must be observed during all phases of operation, service, and repair of this automated staining system. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates the safety standards of the design, manufacturing, and intended use of this system. BioGenex Laboratories, Inc. assumes no liability for the user's failure to comply with these requirements.

I. Precautions and Warning

Xmatrix® Staining System must be unpacked, inspected, and installed by an authorized BioGenex representative. The serial number should be recorded and the accessories checked against the accessory list.



Fuse type: MDA10,250V

Warning: for continued protection against fire hazard, replace only with same type and rating of fuse.



- 1.1 Do not operate the system in an environment where flammable vapors may be present. Operating AC-line-powered equipment in such an environment may cause an explosion.
- 1.2 When replacing fuses on the automated staining system, make sure that the correct fuse is used. Using the wrong fuse may lead to fire or instrument failure.
- 1.3 The AC-line cord provided with the staining system is intended to ground the chassis to help prevent shock and injury to personnel. Xmatrix® meets all national and local electrical codes, which require a three-conductor AC outlet that connects the third wire of the line cord to the earth ground. Do not in any way defeat the grounding connection of this instrument.
- 1.4 Be aware that during or shortly after a high temperature step on the workstation, touching the hot Heater Top(s) may cause burns.




- 1.5 Use of Xmatrix® may require the use of hazardous chemicals. Please refer to the reagent manufacturer's instructions and be sure to follow all applicable regulations for the use, handling, storage, and disposal of any hazardous chemicals.
- 1.6 Appropriate safety precautions should be exercised to protect instrument operators from the risk of biological hazards. At a minimum, gloves should be worn during the operation of the instrument. Gloves, however, may not be sufficient personal protection in all cases. It is the operator's responsibility to confirm with their safety officer what precautions are required.

- 1.7 Liquids used on Xmatrix® System may present a slip hazard if spilled on the floor. Please ensure that all the bulk fluidics and waste containers are placed in a secondary containment.
- 1.8 Liquid waste produced from running the Xmatrix® system may be hazardous and should be disposed off properly according to local regulations.
- 1.9 Excessive jarring of the workstation while the housing is open may cause the housing to close unexpectedly.
- 1.10 Test for ground continuity between the system chassis and adjacent equipment or metal plumbing to ensure that the operator cannot become a conductor between power ground and building ground. Never tamper with the 3-pin power plug. A GFI AC plug receptacle is highly recommended
- 1.11 Do not service the system while power is ON. Refer all servicing to a qualified technician.
- 1.12 Accessing the work platform of the workstation while the robotic arm is in motion may cause bodily injury. Ensure that the run is paused prior to accessing the work platform. The workstation is equipped with a door switch which, when open, will stop the movement of the robotic arm.



- 1.13 This symbol identifies the location of the earth grounding point on the instrument.



- 1.14 This symbol alerts the user that infectious or biological hazardous materials may be present in the instrument and the appropriate precautions must be followed before handling the contents of the instrument or disposing of waste materials.



II. Disclaimer Regarding Staining Results

Xmatrix[®] Staining System may be used only for Research Use only as specified in this manual and for appropriate user-defined purposes.

Any other use is prohibited.

Interpretation of any staining results from the operation of the system is solely the responsibility of the user.

III. Precautions Regarding Biohazardous Materials and Waste

Do NOT dispose of hazardous waste down the sink, sanitary sewer, a rainwater runoff or storm drain. Dispose in accordance with local, state and federal regulations. Should the waste carboy become unusable, it too must be disposed of in the proper manner.

Any business or industrial facility, including research and academic laboratories, that generate hazardous waste are required to comply with EPA hazardous waste regulations that are contained in Title 40 Code of Federal Regulations (CFR) Parts 190-399. It is ultimately the responsibility of the end user of *Xmatrix*[®] to determine whether their waste is considered hazardous or not.

For the purpose of the above disclaimer, biohazardous waste is broadly defined as all biologically contaminated waste that could have the potential to cause harm to humans, domestic or wild animals, or plants. Specific examples of biohazardous waste include human/human tissues, blood, or fluids; cell cultures; or human/animal tissues containing infectious agents or recombinant DNA.

Federal and state laws stipulate that each individual who generates hazardous waste is personally liable and is responsible for assuring compliance with regulations and proper hazardous waste management. Please contact your Biological Safety Management Office to obtain full information on your respective procedures to ensure precautions when handling biohazard materials and all biohazardous waste generated on this instrument.

SECTION 3

ON-SITE INSTALLATION

I. Installation Requirements

Installation Personnel

A BioGenex representative or an authorized technician will install the Xmatrix®. All required assembly and connections will be performed at the time of installation.

Installation Site

The installation site must provide adequate access, ventilation, floor loading and power for the Xmatrix®. These requirements are specified below and in other locations of this manual. If in the future it becomes necessary to move the Xmatrix® to another location the following requirements must be considered:

- a) Locate the Xmatrix® within three (3) feet of a suitable grounded electrical outlet.
- b) The Xmatrix® requires a minimum clearance of 18" on either side and 6" clearance in the rear to allow adequate ventilation and access to the main System ON/OFF switch.

WARNING: This access is important in case of an emergency to allow disconnecting the instrument either by use of the main switch, or by the removal of the power cord either from the power inlet on the instrument or from the power supply on the wall. Do not place anything next to the Xmatrix® that may obstruct access to the main switch or power cord.

NOTE: The Xmatrix® is not routinely turned ON and OFF by using the main power switch. The instrument and computer can be turned ON and OFF independently.

- c) It is not recommended that the Xmatrix® be located on a carpeted floor.
- d) When selecting a location, be aware not to place a load upon any floor of the premises exceeding the floor load per square foot area, which such floor was designed to carry. Point loading must also be considered. The leveling pads of the Xmatrix® can apply a point load of approximately 55 lb/in².

- Xmatrix® System is available as 115V (60Hz, 900 Watts) model and 230V (50Hz) model.
- Make sure to connect your instrument to the appropriate power source when installing or operating it.
- Check the voltage setting on the back of the internal computer inside each Xmatrix®. Staining System to confirm that the setting (115V or 230V) matches the local voltage

II. Installation of Buffer Carboys to the Liquid Panel.

There are 7 carboys each(4L), one for IHC application and the other for FISH application.

The complete set of carboys consists of nine (9) Carboys (Figure 6): Two (2) 15L carboys are used for collecting the waste generated, and seven (7) 4L Carboys for seven (7) different solutions. The colored liquid sockets on this panel are presented in this order:

- Red
- Blue
- Orange
- Yellow
- Green
- White
- Black

The corresponding seven (7) buffer carboys for IHC are labeled as:

- IHC (red)
- SS (blue)
- ISH (orange)
- DEPC (yellow)
- DI-Water (green)
- Alcohol (white)
- Dewax (black)

The corresponding seven (7) buffer carboys for FISH are labeled as:

- Buffer-1 (red)
- Buffer-2 (blue)
- Buffer-3 (orange)
- Buffer-4 (yellow)
- DI-Water (green)
- Alcohol (white)
- DeWax (black)

Before attaching the buffer carboys, make sure the caps are secured tightly to the carboys. Make sure that the carboys caps are not over tightened also. Over tightening may result in cracking of the caps.



Set of 7 Buffer Carboys

The list of carboys for IHC and FISH are listed in the **Table 1a and 1b** respectively. There are seven (7) pairs of sockets on the Liquid Panel. Each pair consists of a white connector for air tubing, and a colored connector for liquid tubing connection (Figure 1a & 1b). When connecting carboy to liquid panel, be sure to match the color of a carboy connector with the color of a connector on this liquid panel. Please connect the carboys as described in **Table 1a and 1b**.

Figure 1a: “Liquid Panel” of the Workstation




Figure 1b: “Liquid Panel” of the Workstation




Fig 03: Loading/Unloading a Carboy

To install Carboys:

Xmatrix® Seven Liter (7L) Carboy Installation Instructions

The 1st step to insert a carboy and lock with locate . Fig.3. press the tab to release the carboy form lock. Insert carboys according to the color labels by sliding carboy into the appropriate position. Push into position until locked into the connectors (Fig.2). Ensure the pressure returns to >5.0 on the pressure display (Fig.4).



Figure 2: Loading Carboys



Figure 4: Pressure Gauge

To release carboys for filling or cleaning, pull the release handle above the carboy (Fig.3). Pressure and lock will release and carboy can be removed (Fig.5).



Figure 5: Remove Carboy

Please connect the carboys as described in Table.1a for IHC run and 1b for FISH run.

Table 1a: Connection Method of 7 Carboys to the Workstation for IHC

Carboy Order in the Shelf (From Left to Right)	Connector Color of Carboy	Type Of Tubing	Socket Color on Manifold
IHC Wash	Red	Black Viton [®]	Red
SS Wash	Blue	Black Viton [®]	Blue
ISH Wash	Orange	Black Viton [®]	Orange
DEPC Water	Yellow	Black Viton [®]	Yellow
DI Water	Green	Black Viton [®]	Green
Alcohol	White	Black Viton [®]	White
DeWax	Black	Black Viton [®]	Black

Table 1b: Connection Method of 7 Carboys to the Workstation for FISH

Carboy Order in the Shelf (From Left to Right)	Connector Color of Carboy	Type Of Tubing	Socket Color on Manifold
Buffer-1	Red	Black Viton [®]	Red
Buffer-2	Blue	Black Viton [®]	Blue
Buffer-3	Orange	Black Viton [®]	Orange
Buffer-4	Yellow	Black Viton [®]	Yellow
DI Water	Green	Black Viton [®]	Green
Alcohol	White	Black Viton [®]	White
DeWax	Black	Black Viton [®]	Black

NOTE: When connecting, insert the 4L Carboy into matching color slot, push until connectors are seated firmly into manifold. Repeat the same for other Carboys. Two (2) waste bottles are placed on the lower deck in the front cabinet of the cart.

CAUTION

Xmatrix® Staining System should be installed by authorized BioGenex service personnel. Once installed and tested on site, the system is intended to remain in that location. If it is necessary to relocate the staining system, contact BioGenex to prepare the system for transport and reinstallation.

Failure to properly prepare the system for transport may result in damage to the automated mechanisms of the system.



This symbol alerts the user that appropriate precautions must be followed when connecting carboys to the cart. Be sure to match the color of each carboy connector with the color of the corresponding socket in the manifold when plugging in the connectors.

Part 2: Priming of Second 2X SSC Carboy During Run Initiation

After initiating a run on the Xmatrix system, the Verifying Consumables and Prime Screens will be displayed. At the “Prime Solutions” screen (**Figure 11**), the Buf. 3-(Orange) box will not be checked. In order to prime the second 2X SSC solution in the carboys, the operator must manually place a check mark in the box next to the “Buf. 3-(Orange)” (**Figure 12**) and continue.



Figure 11: Prime Solutions Screen




Figure 12: Prime Solutions Screen with Buf. 3-Orange checked for priming

III. Waste Carboy Connection

Install the Waste Carboys by attaching the volume sensors and waste tubing. When the waste tubing is connected to the waste carboy, make sure that the connector snaps into the mating connector on the waste carboy. Failure to do so may result in the waste backup in the tubing and may cause flooding of the instrument platform.




Figure 13a: Waste Carboys




Figure 13b: Waste Carboys




Figure 13c: Waste Carboys



Figure 14: Volume Sensor




Figure 15: Waste Tubing

IV. Glossary

Term	Definition/Description
Bulk Bottles	Carboy for containing special liquids such as buffer and DI water. A complete set of buffer carboys for one (1) workstation consists of seven (7) buffer carboys.
Coverslip Mechanism	A combination of the Coverslip (CS) Box and the coverslip tray. The box is the CS container, and the tray holds the boxes in place. There are six (6) CS boxes in three (3) pairs of different sizes (18x18, 25x25, 25x40 mm) in a workstation.
E-Stop Button	Emergency stop button, located at the front of the Xmatrix® workstation.
Gantry	Moving mechanical stand to which the robotic arm is attached and moves back and forth over the working platform.
Housing & Cabinet	Metal hood that encloses the entire Xmatrix® workstation and its accessories. The entire workstation sits on top a cabinet. All buffer carboys and waste carboys fit in the cabinet. The cabinet can be moved around on wheels, which lock into position once the unit position is fixed.
Large Pipette Tip Tray	Blue removable tray with 96 slots for a total of 96 BioGenex supplied blue pipette tips (1 ml each).
Liquid Panel	Connection panel visible on the front side of the cabinet. The Liquid Panel is loaded with seven (7) connection sockets for connecting buffer carboys to the workstation.
Oil Bottle	Oil bottle mounted on the robotic arm.
Oil Pen (Oil Needle)	Slender syringe attached to the lower end of the robotic arm. The Oil Pen applies oil around the slide barrier. The size of the region bounded by oil is determined by the choice of cover slip size.
Pipette Adapter	Pipette adapter on the Z-Head is used to pick up, transport and discard a pipette tip.
Pipette Well	Vertical opening, surrounded by four (4) walls, at the front left part of the work base. A Waste Pipette Tray should be placed under the Pipette Well for receiving the used pipette tip(s) discarded by the Z-Head.
Pressure Gauge	Displays the present air pressure and vacuum pressure in the workstation. The gauge is located at the right end of the Liquid Panel.

Term	Definition/Description
Reagent Tray	A compartment in the front left corner of the workstation for holding reagent vial(s). Up to 49 reagent vials can be accommodated.
Robotic Arm	XYZ robotic arm allowing tri-axial movement with a unique multifunctional staining head that can pick up and discard pipette tips, apply reagent to tissue sections, rinse the slide(s), and blow the excess buffer off slide(s) at the end of a wash step. The Suction Cup, Robotic Head, Oil Pen, and Pipette Adapter are mounted on the robotic arm.
Robotic Head	Function module attached to the Z-Head that sprays liquid onto a slide and/or blow away liquid from a slide.
Slide Block	The largest separated space on the processing platform inside the workstation. A total of four (4) Slide Carriers are placed in this room, and each of the Slide Carriers contains ten (10) Slide Holders. Therefore, it allows up to 40 slides to be processed per run.
Slide Carrier	A long and rectangular compartment that holds ten (10) Slide Holders.
Slide Holder	A rectangular compartment on a Slide Carrier.
Small Pipette Tip Tray	Blue removable tray with 96 slots for a total of 96 BioGenex supplied white pipette tips (200 µl each).
Suction Cup	Translucent rubber cup for picking up Coverslips from a Cover-Slip Holder, or from a slide. The Suction Cup is attached to the lower back end of the Z-Head.
Waste Carboys	Two (2) carboys of 15 liter capacity each are stored in the cabinet at the base of the system for waste reagent from the system.
Waste Coverslip Box	Stainless steel box that collects used Coverslips.
Waste Pipette Tray	Removable tray, placed under the Pipette Well, for holding the used pipette tip(s) discarded by the Z-Head during operation.

SECTION 4**STARTING THE SYSTEM****I. Loading Consumables**

Step	What to Load	Where to Load	How to Load	Precautions
4.1.1	Slide(s) with specimen. Up to <u>40 Slides</u> may be loaded in a single run.	<u>Slide Holders</u> in the 4 Slide Carriers in the Slide Block	a) Lift the latch of the Slide Holder. b) Tilt up the open end of the frame of the Slide Holder with your fingers. c) Insert a Slide into the groove of the Slide Holder. Be sure to place the frosted end of the slide upwards. d) Press the frame into place, and push the Lock button forward to lock slide. e) Repeat the above steps to load other slides into the Slide Carrier(s) as needed. NOTE: Alternatively the slides may be loaded by: 1. Removing the Slide Carrier(s) from Slide Block 2. Loading slide(s) 3. Placing the Slide Carrier(s) back into the Slide Block 4. Turn holder tabs on system to secure Slide Carrier(s)	Make sure the marked end, i.e., the end with the word "BioGenex", of the slide is placed toward the rear of the workstation when you load a slide into the Slide Holder. Ensure slide is loaded sample side up.
4.1.2	Coverslip(s)	Up to <u>6</u> Cover-Slip Boxes at the back left corner of the workstation	1. Load coverslips into coverslip boxes. 2. Attach the coverslip boxes to the coverslip tray in their respective magnetic grooves.	The Cover-Slip(s) should be stacked and loaded in the appropriate tray. Ensure open ends are to the front and back of system, not to the sides.

Step	What to Load	Where to Load	How to Load	Precautions
4.1.3	Reagent(s)	<u>49 Reagents</u> Vials are loaded into the Reagent Tray accordingly	<p>Place vials, one by one into the Reagent Tray without removing the tray from the workstation.</p> <p>NOTE: Always place the top opening of every vial toward the left wall of the Reagent Tray.</p> <p>NOTE: Alternatively the reagent tray can be loaded by:</p> <ol style="list-style-type: none"> 1. Pulling up the Reagent Tray from its seat in the workstation 2. Placing the vials into the Reagent Tray 3. Returning the tray to its seat in the workstation 	Use only BioGenex reagent vials. Incorrect vials could result in fluid detection errors.

Step	What to Load	Where to Load	How to Load	Precautions
4.1.4	Pipette Tips of different volumes (1 ml & 200 µl)	Up to <u>4</u> Pipette Trays of 96 slots in each tray. There are 2 types of Pipette Trays.	<ul style="list-style-type: none"> a) Pull up the blue Large Pipette Tray by grasping its two flaps. Remove it from the workstation. b) Fill in or replace with new box of blue pipette tips (1 ml/each tip). c) Place the Large Pipette Tray back and press it down firmly into its position in the workstation. d) Repeat above three steps to load the second Large Pipette Tray. e) Remove a black metal frame of small pipette tray from workstation. f) Insert a whole box of the white Small Pipette Tray with white pipette tips (200 µl/each tip) in it. g) Place the black metal frame of small pipette tray back into its position in workstation. h) Repeat above three steps to load the second Small Pipette Tray. 	Use only BioGenex supplied pipette tips. Other pipette tips may cause malfunction of the workstation.
4.1.5	Mineral Oil	Mineral Oil Reservoir on robotic head	Call BioGenex Customer Service to fill oil bottle	



Fig: A Coverslips Loading




Fig: B Coverslip stack placing in to respective grooves



Figure C Working Plat Form

II. Starting the System Flowchart






Fig 16: Workstation



Figure 17: Select Username and enter Password in the “Log on” window. There is no password for the Administrator account by default. A password can be set later for this account.




Figure 18a: Selection **Window** offers four options. Select the “Run Xmatrx Application” option.




Figure 18b: Selection **Window** offers two options.
Module A: IHC/ ISH/ SS Protocols, Module B: FISH Protocols

- *Make sure the slides are loaded onto the slide platform.*
- *Make sure the necessary reagent vials are placed on the reagent rack.*

Note: From here onwards Module A refers to IHC, ISH and SS Protocols and Module B refers to FISH Protocols in this manual

SECTION 5

MODULE A – IHC, ISH & SS Protocols

CONFIGURATION SETTINGS

To properly manage a run, and to accurately collect information about the run, a series of configuration settings must be set by the *Administrator* prior to executing the staining run. There are four (4) items within the Configuration module, which will be explained in this section.

I. Maintaining User Profiles

When the *Xmatrix*® system is started for the first time, an authorized *Administrator* will be the only person who can access the software and operate the system. Any other user, who expects to operate the system, must have his/her profile added into the system beforehand by the authorized *Administrator* via the **<User Manager>** window (Figure 20).

1. Select **<Run System Configuration>** from the *Xmatrix* **<Selection Window>**
2. The **<Xmatrix System Configuration & Calibration>** window appears.
3. To access the **<User Manager>** window, click on the **<Configuration>** bar and select **<User Manager>** from the dropdown menu.
4. The **<User Manager>** window will appear.

In the **<User Manager>** window, the *Administrator* can create, modify, or delete a user profile. Related items are User ID, Name, Password, and desired language. The *Administrator* may also limit the protocol type(s) for a user account.

NOTE: The password may also be changed by using the dropdown menu from the **<Edit>** bar in the **<Xmatrix: Staining System>** window.

1. Select the **<Change Password>** option from the **<Edit>** dropdown menu.
2. The **<Change Password>** box (Fig.13a) appears and changes may be made.

The image shows a software dialog box titled "Change Password". It features a blue header bar with the title and a yellow close button. The main area is light blue and contains three text input fields for "Current Password:", "New Password:", and "Confirm Password:". At the bottom, there are two buttons labeled "Save" and "Cancel".

Figure 19 : <Change Password> Window

User Manager ✕

Administrator

Select a User from the List to Modify User Profiles. You can also Delete a User or Create new Users in this screen.

User Profile

* User ID: Administrator * First Name: gg

Language: English * Last Name: gg

User Access

Allow Factory User Protocols
 Allow All Factory and Factory User Protocols
 Allow All User Defined Protocols
 Allow Protocols from List

Figure 20: <User Manager> Window

Xmatrix® permits the selection of different colors, to identify different types of objects in the system, such as Caption Bar and Caption Text, or to identify different status of the same object such as SlideHeat, SlideCool and SlidePause. This rich color scheme helps the user to recognize and distinguish different objects correctly and easily in the course of a run. After entering the required information, click on **<Save Changes>** to save the input to the user database.

II. Specimen Name

Any slide-specific data entered via the **<Specimen Entry Form>** window in the system will be accessible in the **<Slide Editor>** window (Figure 21).

The Slide Editor window is divided into several sections:

- Protocol Information:** Description (IHC PROTOCOL), Number (TMP), Prefix (AX), Kit (IHC Protocols), Group (USER).
- Slide Information:** Specimen Type (dropdown), PreWash (No Wash), Physician Name (dropdown), Patient ID (text), Case Number (text), and a Pattern image.
- Table:** A table with 13 columns: Index, Description, Vol (µl), Cat. Number, Reagent Name, Apply Oil?, Apply CSlip, Heat (°C), Heat (hh:mm:ss), Incubate (hh:mm:ss), Wash Type, Count, Critical Step?, and Pause After?. The table lists steps like BAKING, DEWAX, ANTIGEN RETRIEVAL, etc.
- Buttons:** Expand, Collapse, New Reagent, Print, New Step, Move Up, Copy Step, Save Slide, Cancel, Apply Volume to Reagents (with a dropdown set to 200), Delete Step, Move Down, Paste Step.

Index	Description	Vol (µl)	Cat. Number	Reagent Name	Apply Oil?	Apply CSlip	Heat (°C)	Heat (hh:mm:ss)	Incubate (hh:mm:ss)	Wash Type	Count	Critical Step?	Pause After?
	▶ BAKING												No
	▶ DEWAX												No
	▶ ANTIGEN RETRIEVAL												No
	▶ PEROXIDE BLOCK												No
	▶ ANTIBODY												No
	▶ SUPER ENHANCER												No
	▶ POLYMER HRP												No
	▶ DAB WORKING SOLUTION												No
	▶ HEMATOXYLIN												No
	▶ CLEARMOUNT												No
	▶ XMOUNT												No

Figure 21: **<Slide Editor>** Window

To input specimen names into the system, use the **<Specimen Entry Form>** window (Figure 16a).

To open the **<Specimen Entry Form>** window:

1. From the **<Xmatrix System Configuration & Calibration>** window, select **<Specimen Names>** in the dropdown menu under the **<Configuration>** bar.
2. The **<Specimen Entry Form>** window (Figure 22) will appear.

To input a specimen name:

1. Type in a specimen name into the **<Name>** field. Specimen name field has a limit of 20 characters
2. Select **<Add>**.
3. The name will be listed in the name box below the **<Name>** field of the **<Specimen Entry Form>** window.

4. Select **<Done>** to close the window.

To delete a specimen name:

1. Click to highlight the target name in the name box of the window
2. Select **<Delete>**. The name will disappear.
3. Select **<Done>** to close the window.




Figure 22: <Specimen Entry Form> Window

MODULE B – FISH Protocols

CONFIGURATION SETTINGS

To properly manage a run, and to accurately collect information about the run, a series of configuration settings must be set by the *Administrator* prior to executing the staining run. There are three (3) items within the Configuration module, which will be explained in this section.

Maintaining User Profiles

When the Xmatrix® system is started for the first time, an authorized *Administrator* will be the only person who can access the software and operate the system. Any other user who expects to operate the system, must have his/her profile added into the system beforehand by the authorized *Administrator* via the **<User Manager>** window (Figure 24).

5. Select **<Run System Configuration>** from the Xmatrix **<Selection Window>**
6. The **<Xmatrix System Configuration & Calibration>** window appears.
7. To access the **<User Manager>** window, click on the **<Configuration>** bar and select **<User Manager>** from the dropdown menu.
8. The **<User Manager>** window will appear.

In the **<User Manager>** window, the *Administrator* can create, modify, or delete a user profile. Related items are User ID, Name, Password, and desired language. The *Administrator* may also limit the protocol type(s) for a user account.

NOTE: The password may also be changed by using the dropdown menu from the **<Edit>** bar in the **<Xmatrix: Staining System>** window.

3. Select the **<Change Password>** option from the **<Edit>** dropdown menu.
4. The **<Change Password>** box appears and changes may be made.



Figure 23: **<Change Password>** Window

The screenshot shows the 'User Manager' window with a blue header and a yellow close button. Below the header is a dropdown menu showing 'Administrator'. A bold instruction reads: 'Select a User from the List to Modify User Profiles. You can also Delete a User or Create new Users in this screen.' The main area is divided into two sections: 'User Profile' and 'User Access'. The 'User Profile' section contains fields for '* User ID:' (Administrator), '* First Name:' (gg), 'Language:' (English), and '* Last Name:' (gg). There are buttons for 'Modify ColorScheme' and 'Change Password'. The 'User Access' section has four radio button options: 'Allow Factory User Protocols', 'Allow All Factory and Factory User Protocols', 'Allow All User Defined Protocols', and 'Allow Protocols from List'. An 'Assign Protocols' button is located to the right of these options. At the bottom of the window are three buttons: 'Delete User', 'Save Changes', and 'Cancel Changes'.

Figure 24: <User Manager> Window

Xmatrix® permits the selection of different colors, to identify different types of objects in the system, such as Caption Bar and Caption Text, or to identify different status of the same object such as SlideHeat, SlideCool and SlidePause. This rich color scheme helps the user to recognize and distinguish different objects correctly and easily in the course of a run. After entering the required information, click on **<Save Changes>** to save the input to the user database.

Specimen Name

Any slide-specific data entered via the **<Specimen Entry Form>** window in the system will be accessible in the **<Slide Editor>** window (Figure 25).

The Slide Editor window is divided into several sections:

- Protocol Information:** Description (FISH PROTOCOL), Number (TMP), Prefix (FS), Kit (FISH PROTOCOL), Group (USER).
- Slide Information:** Specimen Type (dropdown), PreWash (No Wash), Physician Name (dropdown), Patient ID (text), Case Number (text), and a Pattern image.
- Table:** A table with columns: Index, Description, Vol. (µl), Cat. Number, Reagent Name, Apply Dil?, Apply CSlip, Heat (°C), Heat (hh:mm:ss), Incubate (hh:mm:ss), Wash Type, Count, Critical Step?, and Pause After?.
- Buttons:** Expand, Collapse, New Reagent, Print, New Step, Move Up, Copy Step, Save Slide, Cancel, Apply Volume to Reagents (with a 200 µl dropdown), Delete Step, Move Down, Paste Step.

Index	Description	Vol. (µl)	Cat. Number	Reagent Name	Apply Dil?	Apply CSlip	Heat (°C)	Heat (hh:mm:ss)	Incubate (hh:mm:ss)	Wash Type	Count	Critical Step?	Pause After?
	Ethanol Wash												No
	Pretreatment												No
	Protease												No
	Aging Operation												No
	Formalin Fixation												No
	Probe Operation												No
	Stringency Wash												No
	DAPI												No

Figure 25: <Slide Editor> Window

To input specimen names into the system, use the **<Specimen Entry Form>** window (Figure 16).

To open the **<Specimen Entry Form>** window:

- From the **<Xmatrx System Configuration & Calibration>** window, select **<Specimen Names>** in the dropdown menu under the **<Configuration>** bar.
- The **<Specimen Entry Form>** window (Figure 26) will appear.

To input a specimen name:

- Type in a specimen name into the **<Name>** field.
- Specimen **<Name>** field has a character limit of 20 characters.
- Select **<Add>**.
- The name will be listed in the name box below the **<Name>** field of the **<Specimen Entry Form>** window.

9. Select **<Done>** to close the window.

To delete a specimen name:

4. Click to highlight the target name in the name box of the window
5. Select **<Delete>**. The name will disappear.
6. Select **<Done>** to close the window.




Figure 26: <Specimen Entry Form> Window

Changing System Clock for Daylight Saving Time (DST)

In August 2005, The United States Congress passed the Energy Policy Act, which changes the start and end dates of daylight savings time. This law goes into effect in 2007. DST will start three (3) week earlier (2:00AM on Sunday March 11) and will end one week later (2:00AM on Sunday November 4). Congress will evaluate the impact on energy consumption in the United States to determine if the new DST will continue beyond 2007.

The Xmatrix® system automatically manages DST settings, but management is based on the DST dates prior to 2007 and controlled by the operating system software. Some releases of system software do not contain the necessary information to autocorrect for the policy changes. Therefore, the automatic management of DST will be disabled, and it will be necessary to manually change the time on the date region that has been set DST to begin and end.

NOTE: Time will need to be changed two (2) times per year until software is upgraded.

Procedure:

1. Double click the time displayed in the Task Bar.
2. Set the correct time from the **<Date & Time>** screen. Press **<Apply>**, followed by **<OK>**.
3. Verify the time is correct.

SECTION 6**MODULE A****BUILDING and INSERTING PROTOCOLS****I. Protocol Types**

There is a protocol template available in the protocol database.

Protocol template:	This protocol exists in the protocol database. This has all the steps included for a typical IHC staining run. User can edit/ modify the protocol parameters. User can also create a new protocol taking this as a template. User can also add or delete the steps in the protocol.
---------------------------	---

II. Staining Steps

On the Xmatrix® system, any staining protocol is a combination of part or all of the following basic staining steps:

Blow	Blows a jet of air over the slide to remove the liquid from the slide.
Incubate	Pauses staining run on the target slide(s) for a specified period of time.
Buffer Slide	Allows buffer to be dispensed onto the particular slide(s) per the applied protocol.
Apply Oil	Applies oil onto the specified slide and seals the four (4) edge lines of a coverslip on the slide.
Apply Coverslip	Delivers a coverslip onto a slide position in the Slide Block.
Remove Coverslip	Removes a coverslip from the target slide and discards the used coverslip.
Heat	Heats up or cools down a slide to the temperature level pre-set in the applied protocol, and maintains the slide at the designated temperature for a specified time.
Repeat	Demands repetition of an existing step or an operation. When this step is present or added, all the other preceding steps listed above within the same operation block in a protocol/slide editor window will be repeated by the specified number of cycles.
Apply Reagent	Uses a 1000 µl pipette tip to apply reagent in volume of 80-850 µl to a target slide.
Apply Micro reagent	Uses a 200 µl pipette tip to apply probe in volume of 10-80 µl to a target slide.

Combinations of staining steps are referred to as *Operations* in the Xmatrix® system. Combinations of different operations may be used to construct a protocol. An Operation can include only one (1) of the following steps:

- Apply Reagent
- Apply Micro reagent

The Xmatrix® System will not allow coverslipping for any operation that contains an “Apply Reagent” step.

III. Building Protocols

Open the “Xmatrix®: **Staining System**” window by selecting the “Run Xmatrix Application” bar in the “Selection Window”




Figure 27: Selection Window

1. Click on the “**Edit**” bar in the “Xmatrix: **Staining System**” window and select **Edit Protocol Database** that leads to **Edit Protocols** window. Click **New** button in this window.
2. In the subsequent window, **New Protocol Editor (Fig.28)** Create a new protocol with the available template. Name the protocol and save it.




Figure 28: New Protocol Editor Window

IV. Inserting a Protocol

Click on a slide position to select it. The selected slide position will be highlighted. To insert a protocol, click on the **Edit** bar in the **Xmatrx Staining System** window and select **Select Protocol** or Right Click the mouse and **Select Protocol**.




Figure 29: Inserting a protocol using the Edit bar




Figure 30: Inserting a protocol using the select 'Select Protocol' command

NOTE: The protocol template currently listed in the software is provided as a reference for use with BioGenex primary antibodies, probes and detection reagents only, and are subject to change. Please consult the appropriate BioGenex product inserts for the most current information before setting up a staining run.

V. Editing Protocols

Only the *Administrator* can edit the default protocol template and other user created protocols, and users can create User Protocols using the <Copy> option.

1. From the <Edit Protocol> window (Figure 31), highlight Protocol template and select the <Edit> button.




Figure 31: <Edit Protocol> Window

2. The <Protocol Editor> screen is displayed (Figure 32). For instructions on using the <Protocol Editor> window, refer to **Section 7: Using Protocol/Slide Editor Windows**.

Figure 32: <Protocol Editor> Window

Copy and Create a new Protocol:

1. From the <Insert Protocol> window (Figure 33), select the desired protocol to copy and select the <Copy> button to open the <New Protocol Editor> window (Figure 34).

Description	Number	Protocol Type	Group
IHC PROTOCOL	AX_TMP	IHC Protocols	USER

Figure 33: <Insert Protocol> <Copy> Window




Figure 34: <New Protocol Editor> Window

2. After completing the design of a new protocol, click on the <Save Protocol> button to save the copied protocol.

Create a New Protocol from an Existing Protocol:

1. From the <Insert Protocol> window (Figure 35), select the desired protocol to copy and select the <New> button to open the <New Protocol Editor> window (Figure 36). For instructions on using the <New Protocol Editor> windows, refer to **Section 7: Using Protocol/Slide Editor Window.**




Figure 35: <Insert Protocol> <New> Window

New Protocol Editor

Protocol Information

Description:

Number:

Prefix:

Kit:

Group:

Slide Information

Specimen Type:

Physician Name:

Patient ID:

Case Number:

PreWash:

Pattern: 

Index	Description	Vol. (µl)	Cat. Number	Reagent Name	Apply Dil?	Apply CSip	Heat (°C)	Heat (hh:mm:ss)	Incubate (hh:mm:ss)	Wash Type	Count	Critical Step?	Pause Alter?
	▶ BAKING												No
	▶ DEWAX												No
	▶ ANTIGEN RETRIVAL												No
	▶ PEROXIDE BLOCK												No
	▶ ANTIBODY												No
	▶ SUPER ENHANCER												No
	▶ POLYMER HRP												No
	▶ DAB WORKING SOLUTION												No
	▶ HEMATOXYLIN												No
	▶ CLEARMOUNT												No
	▶ XMOUNT												No

Figure 36: <New Protocol Editor> Window

2. To save the Protocol, select **<Save Protocol>**. The Protocol can be retrieved later in the **<Insert Protocol>** window.

MODULE B**BUILDING and INSERTING PROTOCOLS**

The user can select any of the options of protocol available for staining runs in the FISH module. The protocol available to the user is,

Default Protocol Template	Default protocol template is given in the software. The user can customize the individual steps/operations in the protocol as per the requirement. This protocol is fully editable and can be created by any user with the permission of creating new protocols.
----------------------------------	--

Staining Steps

On the Xmatrix® system, any staining protocol is a combination of part or all of the following basic staining steps:

Blow Slide	Blows a jet of air over the slide to remove the liquid from the slide.
Incubate Slide	Pauses a staining run on the target slide(s) for a specified period of time.
Wash Slide	Allows buffer to be dispensed onto the particular slide(s) per the applied protocol.
Apply Oil	Applies oil onto the four (4) edge lines of a coverslip on the slide.
Apply Coverslip	Picks up a coverslip from the appropriate coverslip box and applies a coverslip onto a slide.
Remove Coverslip	Removes a coverslip from the target slide and discards the used coverslip into the coverslip dispose box.
Heat Slide	Heats up or cools down a slide to the temperature level specified in the step, and maintains the set temperature on a slide for a time duration specified in the step.
Repeat Step	This step is used in repeating an existing step(s) within the same operation block. In the same operation block, any step(s) listed above this step will be repeated for a number of cycle(s) specified in this step. If there are two repeat steps present within the same operation block, then any step(s) before the first repeat step will be repeated for the number of cycle(s) specified in the first repeat step, and then for the second repeat step, any step(s) between the first and second repeat steps will be repeated for the number of cycle(s) specified in the second repeat step.
Apply Reagent	Uses a 1000 µl pipette tip to apply reagent in volume of 80-850 µl to a target slide.

Apply Probe	Uses a 200 µl pipette tip to apply probe in volume of 10-80 µl to a target slide.
--------------------	---

Combinations of staining steps are referred to as *Operations* in the Xmatrix™ system. Combinations of different operations may be used to construct a protocol. An Operation can include only one (1) of the following steps:

- Apply Reagent
- Apply Probe

The Xmatrix® System will not allow coverslipping for any operation that contains an “Apply Reagent” step.

Create a New Protocol:

1. From the **<Insert Protocol>** window (Figure 37), select the desired protocol and select the **<New>** button to open the **<New Protocol Editor>** window (Figure 38). For instructions on using the **<New Protocol Editor>** windows, refer to **Section 7: Using Protocol/Slide Editor Window**.




Figure 37: <Insert Protocol> <New> Window

New Protocol Editor

Protocol Information

Description:

Number:

Prefix:

Kit:

Group:

Slide Information

Specimen Type:

Physician Name:

Patient ID:

Case Number:

PreWash:

Pattern: 

Index	Description	Vol. (µl)	Cat. Number	Reagent Name	Apply Oil?	Apply CSlip	Heat (°C)	Heat (hh:mm:ss)	Incubate (hh:mm:ss)	Wash Type	Count	Critical Step?	Pause After?
	▶ Ethanol Wash												No
	▶ Pretreatment												No
	▶ Protease												No
	▶ Aging Operation												No
	▶ Formalin Fixation												No
	▶ Probe Operation												No
	▶ Stringency Wash												No
	▶ DAPI												No

Figure 38: <New Protocol Editor> Window

- To save the Protocol, select **<Save Protocol>**. The Protocol can be retrieved later in the **<Insert Protocol>** window.

SECTION 7**MODULE A****USING PROTOCOL/SLIDE EDITOR WINDOWS****I. Editor Windows**

When building or editing a protocol or editing a slide, refer to **Section 6: Building and Inserting Protocols**. There are three types of editor windows available:

New Protocol Editor	Appears when <New> or <Copy> is selected in the <Edit Protocols> window.
Protocol Editor	Appears when <Edit Protocol Database> is selected by clicking the Edit button in the <Xmatrix-Staining System> window.
Slide Editor	Appears by double-clicking on an inserted slide in the <Xmatrix-Staining System> window, or by right-clicking on an inserted slide and choosing the <Edit Slide> option.

The layout of each of these three (3) windows is the same.

Protocol Information	Upper left area of the editor window contains information about the current protocol.
Slide Information	Upper right area which contains information about the current slide (when in the <Slide Editor> window).
Protocol	Central area which contains the operations and steps of the protocol.

The bottom area contains buttons that enable creation or editing of the protocol. Important differences between the three (3) editor windows are:

1. The **<New Protocol Editor>** window allows the user to change the protocol description and number in the **Protocol Information** area, but not the information in the **Slide Information** area. This window allows the user to build a protocol or save a copied protocol.
2. The **<Protocol Editor>** window allows the user to change only the protocol description in the **Protocol Information** area, but none of the information in the **Slide Information** area. This window is used to edit an existing protocol.
3. The **<Slide Editor>** window allows the user to change the information in the **Slide Information** area, but none of the information in the **Protocol Information** area. This window is used to edit single slide only.

The protocol area and the buttons area function identically in all three (3) editor windows.

II. Opening and Closing the <New Protocol Editor> Window

The following information describes how to create a new protocol with the <New Protocol Editor> window. Except where indicated, the techniques used can also be applied to the <Protocol Editor> and the <Slide Editor> windows.

To open the window:

1. From the <Edit Protocol > window (Figure 39), click the <New> button to open the <New Protocol Editor> window (Figure 40).




Fig 39: Edit Protocol Window

Index	Description	Vol. (µl)	Cat. Number	Reagent Name	Apply Oil?	Apply CSlip	Heat (°C)	Heat (hh:mm:ss)	Incubate (hh:mm:ss)	Wash Type	Count	Critical Step?	Pause After?
	▶ BAKING												No
	▶ DEWAX												No
	▶ ANTIGEN RETRIEVAL												No
	▶ PEROXIDE BLOCK												No
	▶ ANTIBODY												No
	▶ SUPER ENHANCER												No
	▶ POLYMER HRP												No
	▶ DAB WORKING SOLUTION												No
	▶ HEMATOXYLIN												No
	▶ CLEARMOUNT												No
	▶ XMOUNT												No

Fig 40: New Protocol Editor Window

2. In the **Protocol Information** area, a protocol number is automatically generated for a new protocol. It can be changed to any three (3) digit number. The description field is blank and must be filled out before the window is closed.

To close the window:

1. After completing the building of a new protocol, click on the button **<Save Protocol>** to save the new protocol and exit the window.

III. Working in the <New Protocol Editor> Window

A protocol is made up of operations, which in turn are made up of steps. An operation is represented by a row with pink cells, and a step is represented by a row with blue cells.

Initially, only the operations are visible.

1. To see the steps of a particular operation, click on the arrow to the left of the operation name.
2. Click the arrow again to collapse the rows.
3. To see all steps for all operations of a protocol, click the **<Expand>** button (Figure 34a) from the buttons area.
4. To hide all steps for all operations, click the **<Collapse>** button.




Figure 41: <New Protocol Editor> window

NOTE: In the <Slide Editor> window, the save button is called <Save Slide> and will save the changes only for that slide.

5. After designing a new protocol:
 - A. Name the Protocol in the description.
 - B. Select <Save Protocol> to save the new protocol.
 - C. The <New Protocol Editor> window will close.
 - D. The <Insert Protocol> window will become active, with the name of the new protocol highlighted.
 - E. By selecting the <Select> button, the <Insert Protocol> window will close. The new protocol will be applied to the Slide Map of the <Xmatrix: Staining System> window.

IV. Button Options in the Editor Windows

Expand	Clicking on this button will display all the steps for all operations on the screen.
Collapse	Clicking on this button will close all the steps for all operations on the screen.

<p>Apply Volume to Reagents</p>	<p>After selecting the volume in the field at the left of this button, click on this button to set the reagent volume in all the <Apply Reagent> steps in the protocol. The reagent volume should not be less than 100 µl. For example, if the selected volume is 95 µl, the selection will not come into effect.</p> <p>An alternative method for specifying the reagent volume in the particular reagent step is to highlight the step and click the cell in the <Vol.> column. User can select the reagent volume from the drop-down menu. The available range for reagent volume selection is 80 to 850 µl. When the <Apply Volume to Reagents> button is clicked, the volume appearing in its field will replace the volume in all the <Apply Reagent> steps.</p>
<p>Apply Volume to Probe</p>	<p>When the selected volume in the related field is 45 or less, the <Apply Volume to Reagents> button will change to <Apply Volume to Probe>. When the <Apply Volume to Probe> button is clicked, the volume appearing in its field will replace the volume in all the <Apply Probe> steps in the protocol.</p> <p>When selecting a number directly in the <Vol.> column, the available volume range for <Probe Step> is 10 to 80 µl.</p>
<p>Print</p>	<p>Clicking on this button will print the protocol displayed in the window.</p>
<p>Save Protocol</p>	<p>Clicking on this button will save a new protocol after the new protocol has been designed.</p>
<p>Cancel</p>	<p>Click this button to cancel all the previous changes made in the window.</p>

V. Slide Information Area

The **Slide Information** area (Figure 35) is open for changes only when working in the **<Slide Editor>** window to edit a single inserted slide. When in the **<New Protocol Editor>** or the **<Protocol Editor>** window, this area is inactive.

When in the **<Slide Editor>** window, the type **<PreWash>** can be selected by clicking on the arrow to activate the dropdown menu.

<Specimen Type> settings can be selected by clicking on the arrow to activate the respective drop-down menus. In order to have these settings available in this window, the specimen types must be entered in advance by an Administrator from the configuration module (see Section 5).




Figure 42: <Slide Information> Area

VI. Protocol Information Area

The <Description> field (Figure 43) is used to enter a title for the new protocol. The <Description> field must be filled in before the Protocol is saved.

The <Number>, <Prefix>, <Kit>, and <Group> fields are automatically filled with relevant information when the window opens. The information in these fields except <Number> cannot be edited by directly working in these fields.




Figure 43: <Protocol Information> Area

VII. Insert Slides After a Run Has Begun through Continuous / Stat Mode

During preliminary login user has an option to use continuous/ stat mode with a check box, As shown in the below Fig CON/STAT log in window. This option need to be selected before starting the run.




FIG: CON/STAT

Once the system has started processing a run, it is possible to insert more slides if needed (and if positions are available on the slide deck) using the <Continuous/ stat> function button.



Note: The Continuous / stat function is available in both Module A / Module B.

slides to the Run through Continuous Mode



Figure 44: <Continuous / stat> button



IMPORTANT: If it is known before the run that this function is likely to be used during the run, the operator should load the reagent carboys to their full capacity, and prime the carboys that are expected to be used prior to starting the run.

To use this function:

1. Select the <Continuous / stat > button at the top of the <Xmatrix: Staining System > window.
2. By Default, Insert Continuous Slides is Checked. Inserting Continuous Slides Option process the slides with Normal Priority .
3. The system will pause to allow insertion of new slides once the current operation has completed. This might take some time depending on the step that is being processed when the <Continuous> function was selected.




Figure 45: Continuous prompt.

4. Read the instructions on the prompt that follows and then select <OK>.

NOTE:



1. The Maximum Number of slides in a staining run should not be more than 20 in Numbers and in three attempts only.

2. Continuous Option Operation may take 30 seconds to Initiate.



CAUTION: After the system pauses, look at the slide map. If any slide is ORANGE (processing) do not proceed! Select the <Cancel> button at the top of the <Xmatrix: Staining System> window, and then <Yes> to confirm the cancellation on the prompt that follows.





Figure 46 : <Continuous / stat> button



NOTE: After cancelling, the operator can retry the <Continuous> option immediately.

5. Once the system has paused, and the operator has verified that no current slides on the slide map are ORANGE, the hood may be opened, and the new slide(s) may be positioned on the slide deck.

NOTE: A <Timeout Occurred> message will appear if <Start> is not selected within 60 seconds.




Figure 47: Timeout Prompt

- 6 Selecting <Yes> on the prompt will allow the operator to continue setting up the new slide(s);
7. However, processing of the existing slides will be delayed until <Start> is selected.
8. Return to the on-screen slide map and right-click on a new slide position to open the <Insert Protocol> window.

NOTE: Multiple slides can be highlighted by using the following Windows options:

Highlight the first slide in a string, hold the <Shift> key, and click the mouse on the last slide in the string •Highlight a slide, hold the <Ctrl> key, and click the mouse on any other slides desired



IMPORTANT: Keep in mind that processing of the existing slides will be delayed until <Start> is selected for the new slide(s).

9. When the desired protocol is highlighted, click <Select>.
10. If the protocol is acceptable as is, proceed to Step 15 below.
- 11 If the protocol needs to be edited, open the <Slide Editor> by double-clicking on a single slide. Refer to Working in the Editor Windows on page 5–25 for instructions on making and saving edits.

12. Verify the slide positions match the slide map and the protocol names and that the correct protocol is chosen for each new slide.

13..Select <Start>. The <Reagent Map> window will appear.

Follow all on-screen instructions and load all required reagents as directed.

VIII. Inserting slides to the Run through Stat Mode.



Figure 48: Continuous / Stat

To use this function:

1. Select the <Continuous / stat > button at the top of the <Xmatrix: Staining System > window.




Figure 49: Continuous / Stat



NOTE: Stat Slides highlighted with Separate color.

NOTE:

1. STAT slides are Processed at the highest priority.
2. STAT Option Operation may take 30 seconds to Initiate.
3. The maximum Number of Slides in a staining Run should be 6 in Number irrespective of the attempts for STAT.
4. Inserted Slides through STAT mode are highlighted with different color.



2. Enable the Insert STAT slides Option and Click <OK>.
3. Insert the slides on to the Available positions in Slide Deck.and click <start>.
4. The system will pause to allow insertion of new slides once the current operation has completed. This might take some time depending on the step that is being processed when the <STAT> function was selected.
5. Once the system has paused, and the operator has verified that no current slides on the slide map are ORANGE, the hood may be opened, and the new slide(s) may be positioned on the slide deck.




Figure 50: Timeout

NOTE: A <Timeout Occurred> message will appear if <Start> is not selected within 60 seconds. (Fig :)

6. Selecting <Yes> on the prompt will allow the operator to continue setting up the new slide(s);
7. However, processing of the existing slides will be delayed until <Start> is selected.
8. Return to the on-screen slide map and right-click on a new slide position to open the <Insert Protocol> window.

NOTE: Multiple slides can be highlighted by using the following Windows options:

Highlight the first slide in a string, hold the <Shift> key, and click the mouse on the last slide in the string •Highlight a slide, hold the <Ctrl> key, and click the mouse on any other slides desired



IMPORTANT: Keep in mind that processing of the existing slides will be delayed until <Start> is selected for the new slide(s).

9. When the desired protocol is highlighted, click <Select>.
10. If the protocol is acceptable as is, proceed to Step 15 below.
11. If the protocol needs to be edited, open the <Slide Editor> by double-clicking on a single slide. Refer to Working in the Editor Windows on page 5–25 for instructions on making and saving edits.
12. Verify the slide positions match the slide map and the protocol names and that the correct protocol is chosen for each new slide.
13. Select <Start>. The <Reagent Map> window will appear.

Follow all on-screen instructions and load all required reagents as directed.

MODULE B**USING PROTOCOL/SLIDE EDITOR WINDOWS****Editor Windows**

When designing or editing a protocol or a slide from the <Insert Protocol> window, refer to **Section 6: Building and Inserting Protocols**. There are three types of editor windows available:

New Protocol Editor	Appears when <New> or <Copy> is selected in the <Edit Protocols> and <Insert Protocols> windows.
Protocol Editor	Appears when <Edit> in the <Insert Protocol> window is selected.
Slide Editor	

The layout of each of these three (3) windows is the same.

Protocol Information	Upper left area of the editor window contains information about the current protocol.
Slide Information	Upper right area which contains information about the current slide (when in the <Slide Editor> window).
Protocol	Central area which contains the operations and steps of the protocol.

The bottom area contains buttons that enable creation or editing of the protocol. Important differences between the three (3) editor windows are:

- The <New Protocol Editor> window allows the user to change the protocol description and number in the **Protocol Information** area, but not the information in the **Slide Information** area. This window allows the user to design a protocol or save a copied protocol.
- The <Protocol Editor> window allows the user to change only the protocol description in the **Protocol Information** area, but none of the information in the **Slide Information** area. This window is used to edit an existing protocol and to modify the staining pattern.
- The <Slide Editor> window allows the user to change the information in the **Slide Information** area, but none of the information in the **Protocol Information** area. This window is used to edit single slide only.

The protocol area and the buttons area function identically in all three (3) editor windows.

Opening and Closing the <New Protocol Editor> Window

The following information describes how to create a new protocol with the <New Protocol Editor> window. Except where indicated, the techniques used can also be applied to the <Protocol Editor> and the <Slide Editor> windows.

To open the window:

- From the <Insert Protocol> window, select the desired protocol and click the <New> button to open the <New Protocol Editor> window (Figure 51).

Index	Description	Vol (µl)	Cat. Number	Reagent Name	Apply Oil?	Apply CSlip	Heat (°C)	Heat (hh:mm:ss)	Incubate (hh:mm:ss)	Wash Type	Count	Critical Step?	Pause After?
▶	Ethanol Wash												No
▶	Pretreatment												No
▶	Protease												No
▶	Aging Operation												No
▶	Formalin Fixation												No
▶	Probe Operation												No
▶	Stringency Wash												No
▶	DAPI												No

Figure 51: <New Protocol Editor> Window

- In the **Protocol Information** area, a protocol number is automatically generated for a new protocol. It can be changed to any three (3) digit number. The description field is blank and must be filled out before the window is closed.

To close the window:

- After completing the design of a new protocol, click on the button <Save Protocol> to save the new protocol and exit the window.

Working in the <New Protocol Editor> Window

A protocol is made up of operations, which in turn are made up of steps. An operation is represented by a row with pink cells, and a step is represented by a row with blue cells.

Initially, only the operations are visible.

- To see the steps of a particular operation, click on the arrow to the left of the operation name.
- Click the arrow again to collapse the rows.
- To see all steps for all operations of a protocol, click the **<Expand>** button (Figure 52) from the buttons area.
- To hide all steps for all operations, click the **<Collapse>** button.

The Slide Editor window is divided into several sections:

- Protocol Information:** Description (FISH PROTOCOL), Number (TMP), Prefix (FS), Kit (FISH PROTOCOL), Group (USER).
- Slide Information:** Specimen Type, Physician Name, Patient ID, Case Number, PreWash (No Wash), and a Pattern selection button.
- Table:** A table with columns: Index, Description, Vol. (µl), Cat. Number, Reagent Name, Apply Dil?, Apply CSlip, Heat (°C), Heat (hh:mm:ss), Incubate (hh:mm:ss), Wash Type, Count, Critical Step?, and Pause After?. The table contains several rows, some expanded to show sub-steps.
- Buttons:** Expand, Collapse, New Reagent, Print, New Step, Move Up, Copy Step, Save Slide, Cancel, Apply Volume to Reagents (with a dropdown set to 200), Delete Step, Move Down, Paste Step.

Index	Description	Vol. (µl)	Cat. Number	Reagent Name	Apply Dil?	Apply CSlip	Heat (°C)	Heat (hh:mm:ss)	Incubate (hh:mm:ss)	Wash Type	Count	Critical Step?	Pause After?
	Ethanol Wash												No
	Wash Slide								00:00:01	Alcohol (White)	1	No	
	Wash Slide								00:02:00	Alcohol (White)	1	No	
	Blow Slide										1	No	
	Pretreatment												No
	Protease												No
	Aging Operation												No
	Formalin Fixation												No
	Probe Operation												No
	Stringency Wash												No
	DAPI												No
	Blow Slide										1	No	
	Apply Micro Reagent	20	111008	DAPI					00:00:01			Yes	
	Apply Coverslip					Yes						Yes	

Figure 52: <New Protocol Editor> Window

NOTE: In the <Slide Editor> window, the save button is called **<Save Slide>** and will save the changes only for that slide.

- After designing a new protocol:
 - Name the Protocol in the description.
 - Select **<Save Protocol>** to save the new protocol.
 - The **<New Protocol Editor>** window will close.
 - The **<Insert Protocol>** window will become active, with the name of the new protocol highlighted.
 - By selecting the **<Select>** button, the **<Insert Protocol>** window will close. The new protocol will be applied to the Slide Map of the **<Xmatrix: Staining System>** window.

Button Options in the Editor Windows

Expand	Clicking on this button will display all the steps for all operations on the screen.
Collapse	Clicking on this button will close all the steps for all operations on the screen.
Apply Volume to Reagents	<p>After selecting the volume in the field to the left of this button, click on this button to set the reagent volume in all the <Apply Reagent> steps in the protocol. The reagent volume should not be less than 100 μL. For example, if the selected volume is 95, the selection will not come into effect.</p> <p>An alternative method for specifying the reagent volume in the particular reagent step is to highlight the step and click the cell in the <Vol.> column. An arrow will display to allow the user to select the reagent volume from the drop-down. The available range for reagent volume selection is 80 to 850 μL. When the <Apply Volume to Reagents> button is clicked, the volume appearing in its field will replace the volume in all the <Apply Reagent> steps.</p>
Apply Volume to Probe	<p>When the selected volume in the related field is 45 or less, the <Apply Volume to Reagents> button will change to <Apply Volume to Probe>. When the <Apply Volume to Probe> button is clicked, the volume appearing in its field will replace the volume in all the <Apply Probe> steps in the protocol.</p> <p>When selecting a number directly in the <Vol.> column, the available volume range for <Probe Step> is 10 to 80 μL.</p>
Print	Clicking on this button will print the protocol displayed in the window.
Save Protocol	Clicking on this button will save a new protocol after the new

	protocol has been designed.
Cancel	Click this button to cancel all the previous changes made in the window.

Slide Information Area

The **Slide Information** area (Figure 53) is open for changes only when working in the **<Slide Editor>** window to edit a single inserted slide. When in the **<New Protocol Editor>** or the **<Protocol Editor>** window, this area is inactive.

When in the **<Slide Editor>** window, the **<PreWash>** and **<Pattern>** settings can be selected by clicking on the arrow to activate the respective dropdown menus.

<Specimen Type> setting can be selected by clicking on the arrow to activate the drop-down menu. In order to have this setting available, the specimen types must be entered in advance by an Administrator from the configuration module.

If desired, the **<Case Number>** may be entered in the appropriate field.





Figure 53: <Slide Information> Area

Protocol Information Area

The **<Description>** field (Figure 54) is used to enter a title for the new protocol. The **<Description>** field must be filled in before the Protocol is saved.

The **<Number>**, **<Prefix>**, **<Kit>**, and **<Group>** fields are automatically filled with relevant information when the window opens. The information in these fields cannot be edited by directly working in these fields.



The screenshot shows a software window titled "Protocol Information". It contains several input fields and a table. The fields are: "Description" (empty), "Number" (007), "Prefix" (FS), "Kit" (FISH PROTOCOL), and "Group" (USER). Below these fields is a table with the following columns: "Index", "Description", "Vol. (µl)", and "Cat. Number". The table contains eight rows of protocol steps, each with a small arrow icon in the "Index" column.

Index	Description	Vol. (µl)	Cat. Number
▶	Ethanol Wash		
▶	Pretreatment		
▶	Protease		
▶	Aging Operation		
▶	Formalin Fixation		
▶	Probe Operation		
▶	Stringency Wash		
▶	DAPI		

Figure 54: <Protocol Information> Area

NOTE:




Inserting the slides to the Run through Continuous / stat function is similar to Module A. In Module B.

SECTION 8

EXECUTING A RUN

I. Understanding a Run



II. Executing a staining run:

The user must ensure that all preparatory steps are completed, such as loading of consumables and reagents. Confirm that all the required buffer Carboys are filled with sufficient volume of appropriate solutions. Ensure the waste carboy is empty. (Refer to **Section 3: On-Site Installation**).

Prior to starting a staining run, the operator should re-verify the following:

1. The slides are properly placed on the slide racks.
2. Required Slide positions are selected.
3. Remove the caps of the reagent vials and load them on the reagent rack with the vial opening on the left. Ensure that the caps are not interchanged when placing them back on the vials.
4. Buffer Carboy levels are sufficient for the run.

Initiating a staining run:

1. Select the slide. Right click/ click twice on the slide to insert the protocol. Click on **<Start>** to start the run.
2. The **<Reagent Screen>** window (Figure 55) will appear. Insert the reagents to be used for the run and click **<Proceed>**.




Figure 55: The Reagent Map Window

3. If required, select **<Print>** in the **<Reagent Map>** window to print the reagent map for future reference.
4. Wait for the workstation to pressurize. The Pressure Gauge at the lower right end of the workstation should display **>5.0**
5. The **<Load Consumables>** window (Figure 56) will appear. Check all consumables and then click **<Proceed>**.

Figure 56: **<Load Consumables>** checklist

Figure 57: **<Select Pipette Positions>** Window

Figure 58: **<Prime Solutions>** Window

6. Make the intended selections in the **<Select Pipette Positions>** window (Figure 57) for the pipette-picking start positions in the two (2) large pipette boxes and two (2) small pipette boxes respectively (front, back or last used position).

7. Select **<Proceed>** to quit this window. The **<Prime Solutions>** box (Figure 58) will appear in order to select the solution(s) that need to be primed.
 8. Buffer lines that are required for the new run will be automatically selected.
 9. Enter the digit(s) into the **<Prime for __ secs>** field to set the prime time (in seconds) (five (5) seconds is the default).
 10. Select **<Start Priming>**. The workstation will start priming.
 11. If priming is not required, select **<Close>** to bypass and the staining run will start.
 12. At the end of priming, the **<Prime Again>** button will become active.
 13. If any of the buffer lines has not been primed sufficiently (no free flow of buffer from the Wash-Blow Head), select **<Prime Again>** again to re-prime solution(s). Or select **<Close>** to proceed to the next step. The **<Prime Solutions>** box will close. The actual staining process will begin. When the process is completed, an alarm will sound.
- 14. Delayed run**
- a. The Delayed Run screen (Fig: 45) will be displayed after clicking on **<close>** on above Screen. If you want to start right away, click **No**.



Figure: 45 Delayed Run

- b. To start a delayed run, click **yes**. The following screen will be displayed.(Fig 46) The system shall ask for time to start processing the slides.



Figure:46 Delayed Run window to set Time and Date

- c. The window will appear to set the time to start the run, and the system shall be able to understand the time set by the user is ahead (later than the current time) and pop up a message to the user when the set time is prior to the current time which is incorrect. Time shall be displayed in 24 Hours format
 - d. You must leave the system on and the program running. Closing the program will cancel the programmed delayed start run even if the program is restarted before the selected start date and time.
 - e. There is a countdown window that appears until the run begins.
 - f. When cancellation of run is selected by the user, system shall ask for a confirmation to cancel to the run.
 - g. Once the set time reaches, system shall start processing the slides without any further response from the user.
15. When the staining run is completed, the **<Run Completed Successfully>** dialog box (Figure 59) will appear for further action.



Figure 59: <Run Completed Successfully> Window

16. Once the run is completed successfully, the user may:
- **<Print Reports>** from the run (refer to **Section 9: Generating Reports** for details of printing)
 - **<Close>** the run
17. If the user selects **<Close>**, a **<Completed>** dialog box appears (Figure 60).
18. Select **<OK>** to close the box and the **<Exit the Application>** bar to exit the subsequent **<Selection Window>**. The operation is now complete.



Figure 60: <Completed> Confirmation Box

III. Monitoring a Run

Real time information of a staining run may be viewed in the <Xmatrix: Staining System> window. The slide map in the window always presents the real time information for each slide being processed. The color scheme together with text on the slide(s), helps the user to easily identify the current step in progress on any slide.

For detailed information on an individual slide:

1. Right click the slide to view the dropdown menu (Figure 61).
2. Select <View Properties> option (Figure 54) to open a <Runtime Properties> window (Figure 62).




Figure 61: <View Properties> Menu



Figure 62: <Runtime Properties> Window for IHC



Figure 63: <Runtime Properties> Window for FISH

The staining progress of a single slide is displayed in the <Runtime Properties> window. Information concerning the current step on a slide as well as the remaining time for that step and slide is presented for quick reference.

NOTE: Estimates of remaining time displayed on the Run time properties window are approximate.

IV. Pausing a Run

To pause a run at any time:

1. Select <Pause> at the top of the <Xmatrx: Staining System> window. The <System Paused> message box (Figure 64) will be displayed.

NOTE: The system will not pause until linked steps are completed.

2. Either click on the **<Resume>** bar in the **<System Paused>** message box, or press the **<Enter>** key on the keyboard to resume the paused staining run.
3. To terminate the run, select **<Cancel the Run>** in the message box.
4. To return the Robotic-Head to its home position, select **<Home Z-Head>**.
5. To park the Z-head, select **<Go To Park Position>**.
6. To remove any completed slides, select **<Remove Completed Slides>**.

CAUTION: Prolonged pause during a staining run may cause invalid test results due to the extended incubation time over the affected staining step(s) or slide(s).





Figure 64: **<System Paused>** Message Box

7. After approximately 60 seconds a **<Confirm>** window will be displayed. If **<Pause>** is not confirmed, the run will resume in 60 seconds.

V. Checking the Availability of Buffers

To check the availability of the buffers during the run:

1. Go to **Operation** in the stainer module



2. Click on Show Buffer Availability, the available volumes are displayed as below:

The 'Buffer Configuration' dialog box displays a table with the following data:

ID	Description	Color	Available Volume(l)
1	XDewax (Black)	Black	0.275000
2	DI Water(Green)	Green	0.450000
3	Alcohol (White)	White	0.225000
4	IHC Wash(Red)	Red	0.550000
5	SS Wash(Blue)	Blue	0.000000
6	ISH Wash(Orange)	Orange	0.000000
7	DEPC Water(Yellow)	Yellow	0.000000

Buttons at the bottom of the dialog include 'Refill Carboy' and 'Done'.

In addition to this, real time monitoring of Buffer level is made available in Bar chart form to the right hand side of <Xmatrx Staining system window>




Figure 52: Real Time Buffer Level Monitoring

In this small Bar Chart, Blue and Red levels show actual carboy level filled. Although buffer levels are calculated and filled as per the requirement, this provides additional information on Live status of Buffer levels. When any level reaches to red color it is indication of buffer has reached to its dead volume level.

VI. Stopping a Run

To stop a run in progress:

1. Select **<Stop>** at the top of the **<Xmatrx: Staining System>** window.

-OR-

2. Click on the **<Operation>** bar to display a dropdown menu. Select **<Exit>** (Figure 65) from the menu.




Figure 65: <Exit> option

3. The stop confirmation message box **<Confirm>** (Figure 66) will be displayed.




Figure 66: <Confirm> Stop Message Box

4. To continue the staining run, click on the **<No>** button. To terminate the current run, click on the **<Yes>** button.
5. Once the **<Yes>** button in the **<Confirm>** box has been clicked, the run will be terminated, and return to the **<Selection Window>** screen.

VII. Emergency Stop

The **<E-Stop>** (Emergency Stop) button is located at the lower right end of the workstation. Pushing this emergency stop button will stop the robotic arm immediately. This is used in case of an emergency. When the **<E-Stop>** button is pressed while the Z-Head is in motion, it stops moving, the system alarm will sound, and a system error message (Figure 67) will be displayed on the screen. Once the **<E-Stop>** button is pressed, the run in progress is cancelled and cannot be resumed.




Figure 67: <System Error Message> Screen

After clearing the emergency situation, release the **<E-Stop>** button that stops the system alarm. Click the **<Close>** button in the system error message window to exit the **<Xmatrix: Staining System>** window. The **<Selection Window>** will be displayed.

When the **<E-Stop>** button is pressed while the Z-Head is not in motion, the system alarm will sound, and the slides on the slide map turn red (Figure 67).




Figure 68: Error Screen

After clearing the emergency situation, release the <E-Stop> button to stop the system alarm. the software exits from the staining application; no need to click on stop button and exit)

VIII. The Typical Error Messages

For most detectable malfunctions in the system, there is a corresponding error message that is displayed on the screen once the malfunction occurs.

The table below lists some typical errors.

The method of clearing an error message is presented along with the error message in the same message window.

No	Error Message	Corrective Action
1	Pipette Boxes for Large Tips are empty	Please Load new boxes and resume the staining run
2	Pipette Boxes for Small Tips are empty	Please Load new boxes and resume the staining run
3	Error picking a small pipette tip from the pipette box	Ensure that tips are available at the desired location and the pipette box is not empty
4	Error picking a large pipette tip from the pipette box	Ensure that tips are available at the desired location and the pipette box is not empty
5	Tip sensor indicates that there is no tip on the Z Head.	Pick up a pipette tip and try again
6	Pipette fell in the middle of an operation	Remove the tip and resume operation
7	Pipette dropped in Slide area	Remove the tip and resume staining
8	Pipette tip has dropped in the reagent rack area	Remove the dropped tip and resume operation
9	Pipette tip has dropped	Remove the fallen pipette tip and resume the operation
10	Unable to dispose pipette	Remove the pipette manually and click <OK> to continue
11	There is not enough reagent available in the reagent vial	Replace the vial or add reagent (enough volume) into the vial.
12	Error sensing the liquid in reagent vial	Make sure enough reagent is available in the vial.
13	The requested pattern is not supported for this operation	Modify the selection of the pattern
14	Large Coverslip boxes are empty	Load coverslips and resume staining on affected slides
15	Error placing the coverslip on the slide	Dispose the coverslip, and place a new coverslip on slide
16	Error disposing the coverslip into the coverslip disposal	Check the calibration and try again. Manually dispose off coverslip.
17	Emergency button is pressed	Click on 'close' to terminate the application.
18	Sensor indicates that the oil bottle is empty.	Please call BioGenex Customer Service.
19	Pipette Drop at reagent vial during aspiration.	Error message pops up and prompts for removing the pipette tip and Resume.
20	Pipette Drop at reagent vial after aspiration	Error message pops up , Prompts for adding Reagent to the vial , click < Yes > to add

No	Error Message	Corrective Action
		reagent if needed. Click < No > . to skip and continue the run.



IMPORTANT: During the pipette tip drop at reagent vial after aspiration of the reagent, Pops up for adding reagent, click < Yes > for adding the reagent to the vial for which the Dispensing of volume is due. Click < NO > to cancel the effected slides.



IMPORTANT: Keep in mind that the Reagent Count down in the Registry takes place when you don't add Reagent to the Run, which may cause insufficient reagent for the run.

SECTION 9

GENERATING REPORTS

I. Launching the Reports Module

There are two (2) ways to access the Reports Module after a staining run has completed:

1. Click the bar **<Run Reports>** in the main dialog window **<Selection Window>** (Figure 69), or
2. Click on the **<Print Reports>** option in the **<Run Completed Successfully>** window (Figure 70).

Either way, the **<Welcome to Reports Module>** window (Figure 71) will appear. This window is the main interface for printing.




Figure 69: <Selection Window> window



Figure 70: <Run Completed successfully> window

II. Default Report Types

The Xmatrx[®] software provides the user with eight (8) types of default report formats. These reports will summarize the staining data for reference. The **<Welcome to Reports Module>** window (Figure 71) displays the various types of reports in the left column. Selecting a report type will display a short description in the box on the right.




Figure 71: **<Welcome to Reports Module>** Window

1. Select the report-type in the **<Select a Report>** dialog. Click on **<Next>** in the **<Welcome to Reports Module>** window, which will open the **<Select Parameters for Run Summary Report>** window (Figure 72).



Figure 72: **<Select parameters for Run Summary Report>** Window

2. Select the date (or date range), run number, or catalog number of the report (Figure 73).

III. Previewing and Printing a Report

Once the selections have been made in the **<Select parameters for Run Summary Report>** window:

1. Select **<Preview & Print Report>** to open the **<Reports>** screen.
2. The **<Reports>** screen displays the report for preview.
3. To print the report, click on the printer icon  at the top of the **<Reports>** screen.




Figure 73: Run Summary Report print view

IV. Exiting the Reports Module

To exit the Reports Module:

1. Click on the **<X>** button on the title bar of the **<Reports>** screen and return to the **<Welcome to Reports Module>** window.
2. Select **<Cancel>** to quit the window, then select **<Exit the Application>** in the **<Selection Window>** to quit Xmatrx®.

SECTION 10

TROUBLESHOOTING

Information in this section lists possible Xmatrix® malfunctions and errors appropriate corrective actions.

More complicated malfunctions or errors may be referred to BioGenex Customer Service.

I. Slide Heating Error

If the system encounters a problem while heating the slide, the process on the specific slide will stop and the slide will be marked as **<Error>** and highlighted in red on the Slide Map screen.




Figure 74: Slide Map: Heating Errors

Corrective Action(s)

1. To view the error message, right click on the failed slide and then click on the **<View Properties>** from the dropdown menu. The **<Runtime Properties>** window displays with the error message displayed on the bottom of the screen in blue.
2. From the Slide Map, resume the failed slide by right clicking on the slide and selecting **<Resume Operation>** from the right click menu. The system will retry the heating step on the affected slide and proceed to the next step if the heating step is successfully completed.

3. If the problem persists, discontinue using the failed slide position and contact BioGenex Customer Service.




Figure 75: Runtime Properties screen

II. Errors at Pipette Tray

If the Z-Head Pipette Adapter fails to pick up a pipette tip from a Pipette Tray when in operation, or if the Pipette Trays are sensed empty, an alarm will sound and an error message **<System Error>** pop-up message will appear for approximately 10 to 15 seconds (Figure 76). Once the **<System Error>** clears, the affected slides are marked with an **<Error>** and highlighted in red in the Slide Map (Figure 77). A **<Fix Errors>** button displays on the left side of the screen. The following steps describe how to recover from this error.




Figure 76: System Error




Figure 77: Fix Errors

Corrective Action(s)

1. Click on **<Fix Errors>** button in the window. The Z-Head will stop proceeding and the **<Fix Errors>** window will appear with the description of the error condition and the appropriate instruction to correct the problem (Figure 78).
 - **<Description>** - Pipette Boxes for Large Tips are empty. Please load new boxes and resume the staining run.
 - **<Fix>** - Refill/replace the empty pipette tip boxes and resume operation.




Figure 78: Fix Errors

2. Perform the action as instructed on the **<Fix Errors>** screen to clear the error condition. In this case, refill the pipette tray, or replace the pipette tip tray with a new tray of pipette

tips. Make sure that the new pipette tray is firmly seated and that there is no obstruction along the path of the Z-Head movement.

3. In the **<Fix Errors>** window, click on the **<Fix Errors>** button.
4. Click on the **<Close>** button to exit from the **<Fix Errors>** window. The system will automatically resume the run as soon as the **<Close>** button is selected.
5. If the problem persists, contact BioGenex Customer Service.

III.A Pipette Drop Error (Before Reagent Aspiration)

If the Z-Head Pipette Adapter drops a pipette tip in the middle of an operation, an alarm will sound and an error message **<System Error>** pop-up will appear on the screen (Figure 79-a)




Figure 79-a: System Error: Pipette Dropped

. Once the **<System Error>** gets addressed by the user, the run will continue.

Corrective Action(s)

Remove the fallen tip and Click on **<OK >** button in the window. When clicked on OK button, then the slides get resumed and the run will continue.

III.B Pipette Drop Error (After Reagent Aspiration)

If the Z-Head Pipette Adapter drops a pipette tip in the middle of an operation after the reagent aspiration, an alarm will sound and an error message **<System Error>** pop-up will appear on the screen (Figure 79-b)




Figure 79-b: System Error: Pipette Dropped

Corrective Action(s)

1. Remove the fallen tip and Click on either **<OK>** or **<Cancel>** button in the window.
 - 1.1 When **<OK>** button is clicked, the lost reagent volume gets restored and the respective slides will be processed. (**Note: prior clicking on Ok button, add the respective volume in the vial manually as asked in the <system error> message displayed**).
 - 1.2 When **<Cancel>** button is clicked, the lost volume doesn't get restored back and the slides which are in error condition remains into error condition without being processed because of the volume insufficiency.

IV. Pipette Disposing Error

After attempting allowable retries, if the Z-Head Pipette Adapter fails to eject a used pipette tip into the Pipette Disposal Tray, an alarm will sound and the error message **<System Error>** will appear on the screen for approximately 10 to 15 seconds (Figure 82). Once the pop-up message clears, the Z-Head will try again to dispose of the used pipette tip into the Pipette Disposal Tray. If the failure recurs (after attempting the allowable retries), an alarm will sound and the error message **<System Error>** will display on the screen for approximately 10 to 15 seconds.




Figure 82: System Error: Pipette Disposing

The system will continue to try and dispose of the pipette tip until the error condition is cleared. If the condition is cleared, the system will normally proceed to the next step. Follow the steps below if the system is unable to clear the error condition on its own.

1. Wait for the next **<System Error>** message pop-up to display. When the **<System Error>** message is counting down and before it expires, manually remove the pipette tip from the Pipette Adapter before the Z-Head starts to move. After the pipette tip is removed and the **<System Error>** is cleared, the system will normally proceed to the next step.
2. If the problem persists, discontinue the run and contact your BioGenex Customer Support.

V. Empty Coverslip Box Error

If the Z-Head suction cup fails to pick up a coverslip from a coverslip stack while in operation, or if the coverslip boxes are sensed as empty, An error message pops with a count down time of 10 sec , the slides in the Slide Map that are to be affected with the error will be marked as **<Error>** and highlighted in red on the Slide Map screen. A **<Fix Error>** button will display on the left side of the screen.




Figure 83: Fix Errors: Empty Coverslip Box

Corrective Action(s)

- Click the **<Fix Errors>** button in the window. The Z-Head will stop proceeding and the **<Fix Errors>** window will appear with the description of the error condition and the appropriate instruction to correct the problem.
 - <Description>** - Coverslip box is empty. Small coverslip boxes are empty. Load coverslips and resume staining on affected slides.
 - <Fix>** - Refill/replace the empty coverslip boxes and resume operation.
- Perform the action as instructed on the **<Fix Errors>** screen to clear the error condition.

Ensure the coverslip boxes are placed in the correct orientation and the coverslips are properly seated in the coverslip box.
- In the **<Fix Errors>** window, click on the **<Fix Error>** button and then click on the **<Close>** button to exit the **<Fix Errors>** window. The system will automatically resume the run as soon as the **<Close>** button is selected.
- If the problem persists, contact BioGenex Customer Service.




Figure 84: Fix Errors: Coverslip Box is Empty

VI. Coverslip Placement / Pickup Error

During the staining operation, if the Z-Head suction cup fails to place a coverslip onto the slide, or fails to pick up a coverslip from the slide, or if it drops a coverslip in the middle of an operation, the slide that is in process will be marked as **<Error>** and highlighted in red on the Slide Map screen.

1. To view the error message, right click on the failed slide and click on the **<View Properties>** from the dropdown menu. The **<Runtime Properties>** window will appear with the error message displayed on the bottom of the screen in blue. See the following two (2) figures for possible error messages.
2. From the Slide Map, resume the failed slide by right clicking on the slide and selecting **<Resume Operation>** from the right click menu. Repeat the step for all other affected slides in the run. The system will retry the coverslip operation on the affected slide and proceed to the next step if it is successfully completed.
3. If the problem persists, contact BioGenex Customer Service.




Figure 85: Coverslip Placement Errors



Figure 86: Coverslip Pickup Errors

VIII. Basic Troubleshooting

No	Symptom	Possible Cause(s)	Corrective Action(s)
1	There is no power to the workstation	<ul style="list-style-type: none"> UPS may not be plugged in Power cord may not be plugged in Power switch not turned ON. No power to the receptacle Safety fuse may be blown 	<ol style="list-style-type: none"> Plug the UPS power cord into an appropriately grounded power source. Plug the power cord into the UPS. Find the power switch on the right side of the workstation, turn it ON. Ensure that the power source has power supply. Request authorized service personnel to check and replace the fuse.
2	Buffers are not primed.	<ul style="list-style-type: none"> Air leakage Air pressure may not be strong enough 	<ol style="list-style-type: none"> Check carboy caps and gaskets, and ensure that they are tightly secured. Also check the connectors on Liquid Panel. Check the Pressure Gauge, and ensure the pressure on the gauge is between 5.2 - 5.7 when priming. If the problem still persists, please call BioGenex Customer Service
3	Failure in Coverslip pick up	<ul style="list-style-type: none"> Related CS Box is empty Vacuum pressure of Suction Cup is not strong enough Suction Cup could not arrive at the correct position over the CS Box 	<ol style="list-style-type: none"> Check the CS Box, and fill in Coverslips as necessary. Ensure there are no leakages in the tubing. Suction Cup could have micro-tear. Call Service to replace the Suction Cup. Request service personnel to calibrate Suction Cup with the CS Box.
4	Failure in pipette tip pick up	<ul style="list-style-type: none"> Target pipette box is empty Pipette Adapter can not fit into a pipette properly 	<ol style="list-style-type: none"> Fill in the Pipette Box, or replace it with a new box full of the correct size pipette tips. Request service personnel to calibrate the Pipette Adapter with the pipette box.
5	Failure in pipette tip disposing	<ul style="list-style-type: none"> Stopping position of Pipette Adapter at Pipette Ejector is not correct 	Request service personnel to calibrate the disposal position of Pipette Adapter at the Pipette Ejector.

SECTION 11

ROUTINE MAINTENANCE

The Xmatrix® System requires minimal routine cleaning and sanitizing. Regular preventive maintenance ensures dependable and consistent staining and helps to eliminate any potential contamination.

I. Daily Cleaning

Daily cleaning should be performed after the last run of the day. It is recommended that the operator wear gloves.

1. Remove and empty the Waste Pipette Tray from the workstation at the end of every workday. Dispose the used pipette tips and rinse the Waste Pipette Tray with DI water. Place the Waste Pipette Tray back into its position in the workstation after cleaning.
2. Use a soft and damp cloth or lint-free paper to wipe out any dirt from all four (4) of the pipette trays (the trays only, not the tips) if necessary.
3. Remove all reagent vials from the Reagent Rack. Remove the Reagent racks and rinse with DI water. Return onto the instrument after drying.
4. Remove and empty the Waste Coverslip Box. Rinse with DI water and dry.
5. Remove all slides from the Slide racks. Remove all four (4) of the slide carriers and soak them in DI water to remove all traces of residue. Dry them with a cloth. Return them to their position.
6. Run System Clean-up (Log in as System Administrator):
 - a. From the <Selection Window>, select <Run System Configuration>.




Figure 87: <Run System Configuration>

- b. The **<Xmatrix – System Configuration & Calibration>** window will appear. Click on the **<Test>** menu.



Figure 88: <Xmatrix – System Configuration & Calibration>

- c) In the **<Test>** menu, select **<System Cleanup>**.



Figure 89: <System Cleanup>

- d) Make sure that the DI water and ethanol carboys contain an appropriate volume to perform the maintenance.
- e) The **<System Cleanup>** window will appear. Click **<Start>**.




Figure 90: <System Cleanup>

NOTE: System Cleanup will complete in approximately 30 minutes. It will utilize the water and ethanol carboys.

- f) After System Cleanup has completed, click on the **<X>** in the upper right hand corner or the screen to close the window.
7. Use a soft and damp cloth to wipe off all of the visible interior and exterior surfaces of the instrument.
- CAUTION:** It is not recommended to wash or rinse any area or surface of the instrument, especially the electrical components and circuits that sit above the Heater Base or under the Inspection Plate of the Housing. This could damage the instrument, causing an electrical short.
8. Use a clean and soft cloth to dry off all the surfaces of the instrument.
 9. Clean the area surrounding the Xmatrix® system.

Prime Solutions

The following are guidelines for operating in the <Prime Solutions> window:

Double-click on **Xmatrix** icon on the computer screen to open <Log on> window

As Administrator, click on <Log on> button in <Log on> window to open <Selection Window>, then click on <Run System Configuration> bar to open the <Xmatrix – System Configuration & Calibration> window

Click <Test> button in <Xmatrix – System Configuration & Calibration> to display a dropdown menu, select the option <Prime Solutions> to open the <Prime Solutions> window.

Operate in <Prime Solutions> window:

- Click in the top four (4) squares to select all (i.e., Dewax, DI water, Alcohol, and Buffer Solutions).
- Click to select one (1) of the four (4) Buffer Lines (repeat this selection with the remaining four (4) Buffer Lines after each round of priming, but deselect the Dewax, DI water, and alcohol after the first pass).
- Click <Start Priming> button to start priming. Repeat the priming twice with a single buffer line. Click <Close> button to end priming with this particular buffer line.
- Repeat above operation to complete the priming with all four (4) buffer lines.
- Connect the alcohol carboy to Dewax line and perform priming. Connect the DI water carboy to all buffer lines except Dewax and Alcohol one after the other and perform priming. Example: Remove the IHC buffer carboy, connect the DI water carboy to the IHC buffer line and perform priming of IHC line. This ensures DI water is run through the IHC line. Repeat this for remaining buffer lines.
- After priming every buffer line with DI water, click <Close> button to end priming.

Click the <X> button to close <Xmatrix – System Configuration & Calibration>, and click <Exit the Application> bar to exit <Selection Window>. The priming operation is done.

II. Monthly Maintenance Schedule

The purpose of monthly maintenance is not only to clean the instrument, but also to prevent the growth of bacteria or mold within the tubing, bulk bottles, and in some areas of the instrument. This procedure should also be followed for disinfecting the instrument prior to transport or storage of the instrument.

1. Perform the daily cleaning procedure.
2. Move the Robotic-Head to the position above the drain hole of the Slide Block. There are three ways to do this:
 - i. Open the **<Prime Solutions>** window in the **<Xmatrix-System Configuration & Calibration>** module (see details in the previous flowchart). The Robotic-Head will move over the drain hole in the Slide Room automatically.
 - ii. Turn OFF the power to the workstation. Slide the gantry forward until the Robotic-Head points to the drain hole of the base plate.

CAUTION: The computer can still be powered ON while the instrument is powered OFF. Never turn OFF the computer while the instrument is ON.

- iii. From **<Xmatrix-System Configuration & Calibration>** select **<Test>** then **<Diagnostics>**. The **<Diagnostics – Motion>** window will open (Figure 91).

In the **<Diagnostics – Motion>** window, use the X/Y arrows to bring the Z-Head at the correct position directly above the drain hole of the instrument.

Or, enter the X 1000 / Y 2100 / Z 0 coordinates in the X, Y, Z fields under **<Current Position>**, then click on the **<Go To Position>** button. The Z-Head will move to the specified position.



Figure 91: **<Diagnostics – Motion>** window

3. Clear off De-Wax residue from De-Wax line: If De-Wax solution has been used through the De-Wax line, empty the De-Wax carboy; fill it with 100 ml of alcohol. Prime De-Wax line three (3) times with the alcohol. Let it sit for ten (10) minutes after the final priming. The alcohol should be rinsed away from De-wax line later. (See below).

NOTE: Skip the above step if the De-Wax feature was not utilized during any staining operations prior to this monthly maintenance schedule.

4. Rinse all the buffer carboys with DI-water: Disconnect all seven (7) buffer Carboy from the instrument. Remove the caps from the buffer carboys and empty each carboy. Fill in each bottle with one gallon of DI-water and add a capful of chlorine bleach. Shake the bleach solution thoroughly in every carboy. Empty the bleach solution from the buffer carboys.

WARNING: Before rinsing the Waste Carboy, turn it in to your Biological Safety Management Office for a proper treatment in accordance with EPA hazardous waste regulations and other local, state, and federal regulations.

5. Prepare fresh cleaning solution for each buffer carboy: Empty the buffer carboys, and use one (1) carboy to dilute 150 ml of Professional Lysol Deodorizing cleaner (the active ingredient: alkyl dimethyl benzyl ammonium chlorides) to 3,600 ml with distilled water.
6. First rinse the liquid transportation system with Lysol Deodorizing cleaning solution: Recap each carboy filled with the above cleaning solution; reattach them to their appropriate ports on the Liquid Panel. Prime their liquid lines three (3) times. Let sit for ten (10) minutes after the final priming. Empty the Carboys.
7. Second rinse of the Carboys with DI-water: Fill each of the seven (7) Carboys with 1 to 2 liters of DI-water only. Re-attach the buffer carboys (i.e., including Waste Carboy) to their appropriate connectors on the instrument. Rinse their liquid lines by priming three (3) times. Empty the Carboys, store them upside down and allow to dry.
8. Clean the drainage tray (i.e., the base of the instrument) with powdered cleanser (i.e., bleach) and a soft plastic brush. Rinse and dry the tray. Ensure that the waste liquid can be drained properly. Check and adjust the drainage tubing as necessary.
9. Dry the interior of the instrument by wiping all surfaces with a clean soft cloth.
10. Leave the Housing Cover open and let the instrument air-dry overnight.

III. Half Yearly Maintenance Arrangement

Only a trained service technician, or BioGenex service personnel, will perform this preventive maintenance every six (6) months. For documentation purposes, it is required that the maintenance report with the instrument serial number, and the maintenance date be clearly filled out. See Appendix 2 for PM checklist

APPENDIX 1

SERVICE PLAN

For a period of one year from the date of installation ("Service Period")*, BioGenex Laboratories, Inc. ("BioGenex") will provide the following services for maintenance and repair of the Xmatrx[®] Automated Staining System ("System") purchased, leased or rented by the Customer ("Services"):

1. Preventative Maintenance Service
2. Emergency Repair Service
3. Software Upgrades

I. Customer Responsibilities

Customer must (1) maintain the installation site and its environment in a condition suitable for operation of the System, and (2) maintain the System under the daily/monthly routine maintenance schedule as provided in the System Operating Manual.

II. Service Charges

BioGenex may charge Customers its standard service rates or may decline to provide the services required to correct a malfunction caused by Customer's failure to fulfill its responsibilities hereunder: the failure of anyone other than BioGenex or its service contractor to comply with its written instructions or recommendations, the combination of the System with an incompatible third party product, the alteration or improper storage, handling, use or maintenance of any part of the System by anyone other than BioGenex or its service contractor, any factor external to the System or beyond BioGenex reasonable control.

III. Additional Terms

Services will be provided during regular business hours. BioGenex may subcontract with service contractors any of the service obligations to Customer. No such subcontract will release BioGenex from its obligations to Customer. In lieu of repair, BioGenex may elect to replace the System or any part thereof, or refund the portion of the Purchase Price paid by Customer. Replacement parts will be provided on an exchange (refurbished) or new part basis, at BioGenex's option. Labor to install replacement parts is included. Replaced parts become BioGenex's property. Components of the System covered by a third party warranty are subject to the limitations contained in that warranty. The Services are due only if the System remains at the Installation Site. Rigging or facility services, accessory, consumable and supply items are charged separately. The Services provided under this Service Plan are performed under the terms and conditions of the Standard System Purchase and License Agreement or Lease and Reagent Purchase Agreement between BioGenex and Customer, as applicable, which are incorporated herein by reference.

IV. Warranty Disclaimer

The System, including all hardware, software, and parts are provided AS IS. BioGenex SPECIFICALLY DISCLAIMS ALL EXPRESS, STATUTORY OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NON-INFRINGEMENT OF THIRD PARTY RIGHTS. Where disclaimer of implied warranties is prohibited by law, such warranties expire 90 days after the date of delivery of the System.

V. Diagnostics and Calibration

BioGenex service personnel or an authorized technician will perform any service, including calibration. An operator should not conduct calibration with this system. When in need, please call BioGenex Customer Service.

The administrator can make a good use of the diagnostics feature of the Xmatrix® system to test the functioning of the system, or to locate the problem with the system. In the Xmatrix® system, there is a diagnostics program ready for use. The most basic module in diagnostic program is the <Test> function. The future **<Diagnostics>**, in which the **<Operation>** screen is an all-inclusive window for testing any of the visible mechanical performances of the Z-Head assembly. Specifically the following operations can be tested:

- Slide operations
- Pipette operations
- General operations
- Reagent operations
- Coverslip operations

To enter **<Operations>** window:

1. Double-click on the **<Xmatrix>** icon on the computer screen to open the **<Log on>** window.
2. Enter as Administrator and enter password, and select **<Log on>** to open the **<Selection Window>**
3. Select **<Run System Configuration>**. The **<Xmatrix – System Configuration & Calibration>** window appears (Figure 93).
4. Select **<Test>** to open a dropdown menu.
5. Select the option **<Diagnostics>** in the dropdown menu to open a submenu.
6. Select **<Operation>** in the submenu. The **<Operations>** screen (Figure 94) will appear.

Upon entering the <Operations> window, five (5) modules will display. By clicking on any button under a module, an operation window will open to work in. Follow the instruction in the window to proceed.




Figure 92: <Xmatrix – System Configuration & Calibration> Window



Figure 93: <Operations> Window

Under the <Diagnostics> dropdown menu, there are tests for <Sensors & Latches>, <Temperature>, and <Motion>. Service personnel will guide the user through these tests, if necessary.

Appendix 2:**Preventive Maintenance Checklist for Xmatrix® System**

S.No.	Items to check	Result (Circle as appropriate)		Note
		Yes	No	
1.	Monitor display is very clear to see	Yes	No	
2.	Keyboard functions well	Yes	No	
3.	Run the "Check Disk C:" feature	Done	No	
4.	Clear previous data of staining runs from computer	Done	No	
5.	Check and update software version of Xmatrix	Done	No	Save data to diskettes if necessary
6.	Check connection of quick connectors and valves on Liquid Panel	Done	No	
7.	Check the overall condition of Z-Head assembly, and tighten the screws	Done	No	
8.	Check Housing Cover Switch and E-Stop switch	Done	No	
9.	Check the connection of all cables and wires to the Electrical Panel	Done	No	
10.	Calibrate Z-Head movement and test Tip Sensor, etc.	Done	No	
11.	Use "Calibration" screen to calibrate Z-Head assembly	Done	No	
12.	Test functionality of Z-Head assembly by Diagnostics	Done	No	
13.	Conduct "Temp Switching Test", "Peak Temp Test", and "Temp Durability Test"	Done	No	
14.	Conduct a "4-Corner Test Run"	Done	No	
15.	Follow monthly maintenance schedule to complete cleaning work	Done	No	Gaskets/O Rings of the caps O.K?
16.	Check the Viton tubing in the carboys and replace if necessary.	Done	No	
17.	Check the gas spring of the lid	Yes	No	

18.	Replace the Assy Carboy 4L Dewax (Part # 65020-30713) with new one	Yes	No	
19.	Check visually for proper oil dispense. Replace oil needle if needed	Yes	No	
20.	Check visually for accurate coverslip positioning	Yes	No	
21.	Check for any shearing or loosening of the Z Head solenoid shafts	Yes	No	
22.	Check for proper operation of the anti-drip valves	Yes	No	No leakage when OFF
23.	Check the oil tubing (PVC) and replace if necessary.			
Comment				

Maintenance Performed by:

Date:

APPENDIX 3

SOFTWARE LICENSE (BioGenex Software)

This software license agreement is a legal agreement between you (either an individual or an entity) and BioGenex Laboratories, Inc. ("BioGenex"), manufacturer of the *Xmatrix*[®] Automated Staining System. This is a license and not an agreement for sale. You obtain no rights other than those granted you under this BioGenex Software License ("License"). By using the *Xmatrix*[®] Automated Staining System on which BioGenex software program(s) have been preinstalled, you are agreeing to be bound by the terms of this license.

I. LICENSE

The term Software means the programs that are installed in and supplied with a BioGenex automated system (System), such as the *Xmatrix*[®] Automated Staining System, and subject to this License, but excluding the Microsoft[®] software which is subject to a separate license.

BioGenex grants you a non-exclusive license to use the Software only in the System in which such Software is initially installed, and only for so long as you are using such System at the Site approved by BioGenex.

BioGenex will provide one copy of the Software with the initial installation of the System, some of which may be pre-installed on the System's hard drive. You may make one backup copy, which is subject to the terms of this License and which must include the copyright notice(s) and any legend(s). All copies of the Software remain the property of BioGenex or its affiliates or suppliers. No title, ownership or other interest in the Software is transferred to you, except as expressly provided herein. You will not modify, reverse engineer, decompile, disassemble, or derive source code from, the Software. In countries which permit reverse engineering notwithstanding the previous sentence, you agree not to conduct such reverse engineering until you have notified BioGenex in writing and given BioGenex a reasonable time to provide such information as you are entitled by law.

Unless BioGenex specifies otherwise, any modifications to, or new releases or versions of, the Software will be provided subject to this License at applicable charges.

II. TRANSFER OF LICENSE

BioGenex will grant a license to the party to whom you transfer possession of the Software (Subsequent Licensee) only when all of the following terms are met:

(a) BioGenex approves such transfer in writing in its sole discretion; (b) the Software, including all copies, and the System are concurrently transferred to the Subsequent Licensee; and (c) the Subsequent Licensee agrees in writing to be bound by the terms of this Software License with respect to the Software transferred. Transfer of the Software terminates your license from BioGenex.

III. INTELLECTUAL PROPERTY

The Software is the confidential trade secret of BioGenex. You agree to take precautions to prevent any disclosure or use of Software (except as specifically permitted herein) consistent with precautions used to protect your own confidential trade secrets, but in no event less than reasonable care. You shall not copy, disassemble, reverse engineer, or decompile the Software.

IV. COPYRIGHT

The Software is owned by BioGenex and is protected by United States copyright laws and international treaty provisions and all other applicable national laws. Therefore, you must treat the Software like any other copyrighted material.

V. WARRANTY

If the Software is warranted, the terms and conditions of such warranty is specified in the agreement under which such Software was originally acquired from BioGenex or a BioGenex Authorized Distributor (Distributor). Otherwise, the Software is provided AS IS. The operation of the Software may not be uninterrupted or error free.

SUCH WARRANTIES REPLACE ALL OTHER WARRANTIES, EXPRESS, STATUTORY OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT OF THIRD PARTY RIGHTS.

Support services for Software which is not under warranty may be available from BioGenex or a Distributor for a charge.

VI. LIMITATION OF LIABILITY

THE MAXIMUM AGGREGATE LIABILITY OF BIOGENEX AND ITS SUPPLIERS, DEVELOPERS, SUBCONTRACTORS, AND AGENTS FOR ANY CLAIMS ARISING OUT OF OR RELATED TO THIS AGREEMENT, WHETHER IN CONTRACT, TORT, PRODUCT LIABILITY OR OTHERWISE, SHALL NOT EXCEED THREE THOUSAND DOLLARS.

Under no circumstances will BioGenex be liable for losses or damages resulting from (1) third party claims against you or claims by you based on third party claims (other than those referred to in item 1 in the first paragraph of this Section), (2) loss of stored, transmitted or recorded data, (3) consequential damages (including lost profits or savings), or incidental, indirect or punitive damages, even if BioGenex or a Distributor is informed of their possibility, or (4) inability to perform under this License, in whole or in part, caused by acts of God or any other event beyond BioGenex's control.

This Section sets forth the maximum collective responsibility of BioGenex, its suppliers, developers, subcontractors and agents, and all such parties are intended beneficiaries of this Section.

VII. U.S. GOVERNMENT RESTRICTED RIGHTS

If you are a U.S. Government user, the Software and documentation are commercial computer software developed exclusively at private expense and are subject to the terms of this License pursuant to 48 C.F.R. 12.212 or 48 C.F.R. 227.7202. The manufacturer is BioGenex Laboratories, Inc., 48810 Kato Road, Suite 200E, Fremont, CA 94538.

VIII. GENERAL

You may terminate your license at any time by destroying all your copies of the Software or as otherwise described in this License. BioGenex may terminate your license if you do not comply with the terms of this License. Upon such termination, you will destroy all your copies of the Software. Any attempt to sublicense, rent, lease, assign or, except as expressly provided herein, to transfer any copy of the Software is void. This License is governed by the laws of the State of California without regard to its rules on conflicts of laws.

APPENDIX 4**IMPORTANT - READ CAREFULLY BEFORE USING THE Xmatrix® AUTOMATED STAINING SYSTEM WHICH CONTAINS MICROSOFT® SOFTWARE****SOFTWARE LICENSE AGREEMENT
(Microsoft Software)**

IMPORTANT-READ CAREFULLY: This End-User License Agreement (EULA) is a legal agreement between you (either an individual or a single entity) and the manufacturer (Manufacturer) of the special purpose computing device (SYSTEM) you acquired which includes certain Microsoft software product(s) installed on the SYSTEM (SOFTWARE PRODUCT or SOFTWARE). The SOFTWARE includes computer software, the associated media, any printed materials, and any online or electronic documentation. By installing, copying or otherwise using the SOFTWARE, you agree to be bound by the terms of this EULA. If you do not agree to the terms of this EULA, Manufacturer and Microsoft Corporation (Microsoft) are unwilling to license the SOFTWARE to you. In such event, you may not use or copy the SOFTWARE, and you should promptly contact Manufacturer for instructions on return of the unused product(s) for a refund.

The SOFTWARE is protected by copyright laws and international copyright treaties, as well as other intellectual property laws and treaties. The SOFTWARE is licensed, not sold.

I. GRANT OF LICENSE

This EULA grants you the following rights:

Software. You may use the SOFTWARE as installed on the SYSTEM.

Application Sharing. The SOFTWARE may contain technology that enables applications to be shared between two or more SYSTEMS, even if an application is installed on only one of the SYSTEMS. You should consult your application license agreement or contact the application licensor to determine whether sharing the application is permitted by its licensor.

Storage/Network Use. If the SOFTWARE PRODUCT is installed on the SYSTEM over an internal network from a server, you must acquire and dedicate a license for the SOFTWARE PRODUCT for each SYSTEM on which the SOFTWARE PRODUCT is used or to which it is distributed. A license for the SOFTWARE PRODUCT may not be shared or used concurrently on different SYSTEMS.

Back-up Copy. If Manufacturer has not included a copy of the SOFTWARE on separate media with the SYSTEM, you may make a single copy of the SOFTWARE for use solely for archival purposes with the SYSTEM.

II. DESCRIPTION OF OTHER RIGHTS & LIMITATIONS

Limitations on Reverse Engineering, Decompilation and Disassembly. You may not reverse engineer, decompile, or disassemble the SOFTWARE, except and only to the extent that such activity is expressly permitted by applicable law notwithstanding this limitation.

Single SYSTEM. The SOFTWARE is licensed with the SYSTEM as a single integrated product. The SOFTWARE may only be used with the SYSTEM.

Rental. You may not rent or lease the SOFTWARE.

Software Transfer. You may permanently transfer all of your rights under this EULA only as part of a sale or transfer of the SYSTEM, provided you retain no copies, you transfer all of the SOFTWARE (including all component parts, the media, any upgrades or backup copies, and this EULA, and if applicable, the Certificate(s) of Authenticity), and the recipient agrees to the terms of this EULA. If the SOFTWARE is an upgrade, any transfer must include all prior versions of the SOFTWARE.

Termination. Without prejudice to any other rights, Manufacturer or Microsoft may terminate this EULA if you fail to comply with the terms and Conditions of this EULA. In such event, you must destroy all copies of the SOFTWARE and all of its component parts.

Single EULA. The package for the SOFTWARE may contain multiple versions of this EULA, such as multiple translations and/or multiple media versions (e.g., in the user documentation and in the software). In this case, you are only licensed to use one (1) copy of the SOFTWARE PRODUCT.

III. UPGRADES AND RECOVERY MEDIA

If the SOFTWARE is provided by the Manufacturer on media separate from the SYSTEM and is labeled "For Upgrade Purposes Only" (Upgrade SOFTWARE), you may install one copy of the Upgrade SOFTWARE onto the SYSTEM as a replacement copy for the SOFTWARE originally installed on the SYSTEM and use it in accordance with Section 1 of this EULA. You may also install additional copies of the Upgrade SOFTWARE as replacement copies onto additional SYSTEMS which are the same brand and model as the SYSTEM and contain a duly licensed copy of the same version and language release of the SOFTWARE (ADDITIONAL SYSTEMS), provided that (1) Manufacturer has supplied a corresponding serialized sticker for each additional copy of the Upgrade SOFTWARE; and (2) you affix a serialized sticker per Manufacturer's instructions for each unit of Upgrade SOFTWARE you install.

If the SOFTWARE is provided by Manufacturer on separate media and labeled as Recovery Media, you may not make a copy of the SOFTWARE as described in Section 1 for archival purposes. Instead, you may use the Recovery Media solely to restore or reinstall the same version and language release of the SOFTWARE as originally installed on the SYSTEM and thereafter use the SOFTWARE as restored or reinstalled in accordance with Section 1 of this EULA. A single unit of Recovery Media may be used by you to restore or reinstall the SOFTWARE on ADDITIONAL SYSTEMS.

IV. COPYRIGHT

All title and copyrights in and to the SOFTWARE (including but not limited to any images, photographs, animations, video, audio, music, text and applets, incorporated into the SOFTWARE), the accompanying printed materials, and any copies of the SOFTWARE, are owned by Microsoft or its suppliers. You may not copy the printed materials accompanying the SOFTWARE. All rights not specifically granted under this EULA are reserved by Microsoft.

V. PRODUCT SUPPORT

Product support for the SOFTWARE is not provided by Microsoft or its subsidiaries. For product support please refer to Manufacturer's support number provided in the documentation for the

SYSTEM. Should you have any questions concerning this EULA, or if you desire to contact Manufacturer for any other reason, please refer to the address provided in the documentation for the SYSTEM.

VI. LIMITED WARRANTY

Manufacturer warrants that the SOFTWARE will perform substantially in accordance with the accompanying written materials for a period of ninety (90) days from the date of receipt. Any implied warranties on the SOFTWARE are limited to ninety (90) days. Some states/jurisdictions do not allow limitations on duration of an implied warranty, so the above limitation may not apply to you.

Customer Remedies. Manufacturer's and its suppliers' entire liability and your exclusive remedy shall be, at Manufacturer's option, either (a) return of the price paid, or (b) repair or replacement of the SOFTWARE that does not meet the above Limited Warranty and which is returned to Manufacturer with a copy of your receipt. This Limited Warranty is void if failure of the SOFTWARE has resulted from accident, abuse, or misapplication. Any replacement SOFTWARE will be warranted for the remainder of the original warranty period or thirty (30) days, whichever is longer.

No Other Warranties. EXCEPT AS EXPRESSLY PROVIDED IN THE LIMITED WARRANTY SECTION ABOVE, THE SOFTWARE IS PROVIDED TO THE END USER AS IS WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF NON-INFRINGEMENT, MERCHANTABILITY, AND/OR FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK OF THE QUALITY AND PERFORMANCE OF THE SOFTWARE IS WITH YOU.

No Liability for Consequential Damages. MANUFACTURER'S SUPPLIERS SHALL NOT BE HELD TO ANY LIABILITY FOR ANY DAMAGES SUFFERED OR INCURRED BY THE END USER (INCLUDING, BUT NOT LIMITED TO, GENERAL, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES INCLUDING DAMAGES FOR LOSS OF BUSINESS PROFITS, BUSINESS INTERRUPTION, LOSS OF BUSINESS INFORMATION AND THE LIKE), ARISING FROM OR IN CONNECTION WITH THE DELIVERY, USE OR PERFORMANCE OF THE SOFTWARE.

VII. NOTE ON JAVA SUPPORT

THE SOFTWARE PRODUCT MAY CONTAIN SUPPORT FOR PROGRAMS WRITTEN IN JAVA. JAVA TECHNOLOGY IS NOT FAULT TOLERANT AND IS NOT DESIGNED, MANUFACTURED, OR INTENDED FOR USE OR RESALE AS ON-LINE CONTROL EQUIPMENT IN HAZARDOUS ENVIRONMENTS REQUIRING FAIL-SAFE PERFORMANCE, SUCH AS IN THE OPERATION OF NUCLEAR FACILITIES, AIRCRAFT NAVIGATION OR COMMUNICATION SYSTEMS, AIR TRAFFIC CONTROL, DIRECT LIFE SUPPORT MACHINES, OR WEAPONS SYSTEMS, IN WHICH THE FAILURE OF JAVA TECHNOLOGY COULD LEAD DIRECTLY TO DEATH, PERSONAL INJURY, OR SEVERE PHYSICAL OR ENVIRONMENTAL DAMAGE.

VIII. U.S. GOVERNMENT RESTRICTED RIGHTS

The SOFTWARE and documentation are provided with RESTRICTED RIGHTS. Use, duplication or disclosure by the Government is subject to restrictions as set forth in subparagraph (c)(1)(ii) of the Rights in Technical Data and Computer Software clause at DFARS 252.227-7013 or subparagraphs (c)(1) and (2) of the Commercial Computer Software -

Restricted Rights at 48 CFR 52.227-19, as applicable. Manufacturer is Microsoft Corporation, One Microsoft Way, Redmond, WA 98052-6399.

If you acquired this EULA in the United States, this EULA is governed by the laws of the State of Washington. If you acquired this EULA in Canada, this EULA is governed by the laws of the Province of Ontario, Canada. If this EULA was acquired outside the United States, then local law(s) may apply.

Should you have any questions concerning this EULA, please contact BioGenex at Tel. No. (800) 421-4149 [U.S. Customers ONLY] or your local BioGenex Distributor.