Doc. No. 932-HK100E-4 Rev. No. J
Release Date: 18-Jul-2024

ENGLISH

Catalog No.	Counter Stain
HK100-5K, HK100-9K	Hematoxylin

S26 S37 S24/25 S60 = In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable gloves. Avoid contact with skin and eyes. This material and its container must be disposed of as hazardous waste.

The following table provides the recommended treatments with counterstain. Refer to the antibody and detection kit datasheets for complete staining protocols. User should validate any other conditions. For automated use follow instructions in operator's manual.

Counter stain	BioGenex Recommendations
----------------------	---------------------------------

Hematoxylin	Stain slides for 1 to 5 minutes, rinse with distilled water. Dehydrate sections in a series of graded alcohols and clear in xylene or xylene substitute, then mount. Compatible with AEC, DAB and Fast Red chromogens.
-------------	--

1. Run appropriate positive and negative controls to ensure the optimal signal.
2. Refer to the troubleshooting section in the package inserts of BioGenex Super Sensitive Detection Systems (or other equivalent detection systems) for remedial actions on detection system related issues, or contact BioGenex Technical Service Department at 1 800-421-4149 or customer.service@biogenex.com to report unusual staining.

Expected Results	
Counter stain	Expected Results
Hematoxylin	Nucleus stains bluish purple.
Methyl Green	Nucleus stains light green
Nuclear Fast Red	Nucleus stains red, cytoplasm stains pink.

BioGenex has conducted studies to evaluate the performance of all its counterstains using several BioGenex antibodies and detection systems. The counterstains have been found to show color contrasts with respect to the primary stain thereby enhancing the staining pattern. BioGenex counter stains have shown reproducible and consistent results when used within a single run, between runs, between lots and wherever applicable between manual and automated runs. The products have been determined to be stable for the periods specified on the labels either by standard real time or accelerated testing methods. BioGenex ensures product quality through 100% quality control for all products released and through surveillance programs.

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Bibliography

- Sheehan DC and Hrapchak BB. Theory and Practice of Histotechnology. St. Louis: C.V. Mosby Co. 1980.
- Kieman JA. Histological and Histochemical Methods: Theory and Practice. New York: Pergamon Press 1981.
- Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and its pitfalls. Lab Med 1983; 14:767-770.
- Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen. A possible source of error in immuno histochemistry. Am J ClinPathol 1980 May;73(5):626-632.

EC|REP
Representative in the European Community
MandatariorellaComunitàEuropea
Bevollmächtigter in der EuropäischenGemeinschaft
Representanteautorizadoen la ComunidadEuropea

IVD
In Vitro Diagnostic Medical Device
Dispositivo medico-diagnostico in vitro
In Vitro Diagnostikum
Productosanitario para diagnóstico in vitro

Consult Instructions for use
Consultarele istruzioni per l'uso
Gebrauchsanweisungbeachten
Consulte las instruccionesde uso

Temperature Limitation
Limiti di temperatura
ZulässigerTemperaturbereich
Limite de temperatura

Manufacturer
Fabbricante
Hersteller
Fabricante