Accelerating the pace of precision medicine

48810 Kato Road, Suite 100E & 200E,

EC REP

Emergo Europe

Westervoortsedijk 60, 6827

AT Arnhem, The Netherlands

Tel : +1 (800) 421-4149,
Fax: +1 (510) 824-1490,

support@biogenex.com

DILUENTS

Doc. No. 932-HK156E-4 Rev. No. I
Release Date: 18-Jul-2024

ENGLISH

Intended Use
BioGenex diluents are intended for in vitro diagnostic use for dilutions, adjustment and optimization of antibody titers or concentrations of reagents for use in an immunohistochemistry (IHC) or in situ hybridization (ISH) assay.

Summary and Explanation
In order to achieve desired staining pattern and intensity the titers of antibodies or concentrations of reagents may need to be optimized. BioGenex offers several diluents for specific antibody or reagent products.

Principles of the Procedure
The diluents have been developed and optimized for use with BioGenex antibodies and reagent products. Optimal dilutions result in enhancement of signal to noise ratio. The diluents help in maintaining the antibody specificity and stability and also prevent contamination.

Reagents Provided		
Catalog No.	Diluents	
HK156-5KE HK156-25K	Common Antibody Diluent	
HK157-5KE	Streptavidin–Peroxidase Diluent	
HK165-5KE	Link Diluent	
HK941-YAKE HK941-25K HK941-4K	Enhanced Antibody Diluent	

Reagents Required but Not Supplied
All reagents required for dilution are provided. See antibody and detection kit datasheets for complete set of reagents required for immunohistochemistry (IHC) or in situ hybridization (ISH) procedures.

Storage and Handling
Store all reagents at 2-8°C. Do not use after expiration dates indicated on the reagent labels.

Precautions
All reagents contain sodium azide at concentrations of less than 0.1%. Sodium azide is not classified as a hazardous chemical at the concentration of this product. However, toxicity information regarding sodium azide at the product’s concentration has not been thoroughly investigated. For more information, a Material Safety Data Sheet (MSDS) for sodium azide in pure form is available upon request.

Dilution of Concentrated Primary Antibodies
BioGenex Enhanced Antibody Diluent is a ready-to-use solution. Dilutions are given as the ratio of concentrated volume to final volume (1:100 means 1 part of concentrated antibody mixed with 99 parts of diluent). For dilutions greater than 1:100, the use of serial dilutions is recommended to improve accuracy. For example, a 1:500 dilution can be achieved by first making a 1:100 dilution of the concentrated antibody and then making a further 1:5 dilution of the 1:100 dilution.

Insufficient dilution of primary antibody (i.e., antibody concentration too high) may cause false positive staining, while over-dilution (antibody concentration too low) may cause false negative staining. Different dilutions of a concentrated antibody is needed starting from the titer recommended by vendors.

Dilution of Concentrated Links and Labels

1. Concentrated link (secondary) antibodies and concentrated labels are shipped with an instruction manual listing the recommended dilutions depending on the level of sensitivity required.
Be sure to match the appropriate diluent with the concentrated label being employed. Peroxidase labels can be diluted in Streptavidin Peroxidase Diluent (DO NOT use Common Antibody Diluent for peroxidase labels since it contains sodium azide, which inhibits peroxidase activity). Alkaline phosphatase labels can be diluted in Common Antibody Diluent.

Staining Protocol
Refer to the antibody and detection kit datasheets for complete staining protocols and recommended dilution factors. User should optimize and validate any other conditions.

Quality Control
Refer to the appropriate detection system package inserts for guidance on general quality control procedures.

Trouble shooting
Refer to the troubleshooting section in the package inserts of BioGenex Super Sensitive Detection Systems (or other equivalent detection systems) for remedial actions on detection system related issues, or contact BioGenex Technical Service Department at (925) 275-0550 to report unusual staining.

Expected Results
Optimal dilution of concentrated antibodies and or reagents should result in deposition of colored chromogen pigment at the site of specific antigen-antibody reaction with minimal to no non-specific background.

Limitations of the procedure
Excessive dilutions may result in false negative or weaker staining. Lower dilution may result in non- specific background staining.

Performance Characteristics
BioGenex has conducted studies to evaluate the performance of all its diluents using several BioGenex antibodies and detection systems. The appropriate dilutions of concentrated antibodies or reagents have been found to show clear pattern and intensity of specific antigen-antibody staining reaction with minimal to no non-specific background. BioGenex diluents have shown reproducible and consistent results. The products have been determined to be stable for the periods specified on the labels either by standard real time or accelerated testing methods. BioGenex ensures product quality through 100% quality control for all products released and through surveillance programs.

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Reagenti forniti	
N. di catalogo	Diluenti
HK156-5KE HK156-25K	Common Antibody Diluent
HK157-5KE HK165-5KE	Streptavidin–Peroxidase Diluent Link Diluent
HK941-YAK HK941-25K HK941-4K	Enhanced Antibody Diluent

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Bibliography

- Sheehan DC and Hrapchak BB. Theory and Practice of Histotechnology. St. Louis: C.V. Mosby Co. 1980.
- Kieman JA. Histological and Histochemical Methods: Theory and Practice. New York: Pergamon Press 1981.
- Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and its pitfalls. Lab Med 1983; 14:767-770.
- Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen. A possible source of error in immuno histochemistry. Am J Clin Pathol 1980 May;73(5):626-632.
- National Institute for Occupational Safety and Health,(NIOSH), Rockville, MD. Explosive azide hazard, Publication No. 78-127, 1976.
- Center for Disease Control. Decontamination of Laboratory Sink Drains to Remove Azide Salts. Center for Disease Control Manual Guide--Safety Management, No. CDC-22, Atlanta, Georgia. April 30, 1976.

	Representative in the European Community Mandatario nella Comunità Europea Bevollmächtigter in der Europäischen Gemeinschaft Representante autorizado en la Comunidad Europea
	In Vitro Diagnostic Medical Device Dispositivo medico-diagnostico in vitro In Vitro Diagnostikum Producto sanitario para diagnóstico in vitro
	Cons ult Instructions for use Consultare le istruzioni per l'uso Gebrauchsanweisung beachten Consulte las instrucciones de uso
	Temperature Limitation Limiti di temperatura Zulässiger Temperaturbereich Limite de temperatura
	Manufacturer Fabbricante Hersteller Fabricante