

**New and Improved Super Sensitive™ Polymer-HRP IHC Detection
for Manual and Automated Systems**

Kit Components:

Cat. No.	No. of Slides	EZ-AR™ 1	EZ-AR™ 2	Peroxid eBlock	Power Block	Mouse Negative Control	Rabbit Negative Control	Super Enhancer™	Poly-HRP	Substrate	Chromogen	Empty Mixing Vial	Hematoxylin	Xmount	Symbol
QD400-60KEN Manual	60	---	---	HK111-5KTN (6 ml)	HK083-5KTN (6 ml)	HK119-5MN (3 ml)	HK408-5RTN (3ml)	HK518-06KN (6 ml)	HK519-06KN (6 ml)	HK520-10KN (10 ml)	HK124-1KN (1 ml)	Mixing Vial (6ml)	HK100-5KTN (6 ml)	---	
QD420-YIKEN Manual	500	---	---	---	---	---	---	HK518-50KN (50 ml)	HK519-50KN (50 ml)	HK520-50KN (50 ml)	HK124-025KN (1 x 2.5 ml)	Mixing Vial (6ml)	---	---	
QD430-XAKEN Manual	1000	---	---	---	---	---	---	HK518-YAKN (100 ml)	HK519-YAKN (100 ml)	HK520-YAKN (100 ml)	HK124-05KN (5 ml)	Mixing Vial (6ml)	---	---	
QD440-XAKEN Manual	1000	---	---	---	---	---	---	HK518-YAKN (100 ml)	HK519-YAKN (100 ml)	---	---	---	---	---	
QD550-YCXEN Xmatrix Infinity	200	HX031-YCXN (16 ml)	HX032-YCXN (16ml)	HX026-YAXN (3x16 ml)	HX083-YAXN (3x 16 ml)	---	---	HX027-YCXN (2x16 ml)	HX028-YCXN (2x16 ml)	HX029-50XN (4 x 13 ml)	HX010-YCKN (4ml)	HX616-YCXN (4no)	HX030-YAXN (3 x 16 ml)	HX035-YCDN (15ml)	
QD550-YCDEN Xmatrix Elite	200	HX031-YCDEN (16 ml)	HX032-YCDEN (16 ml)	HX026-YADEN (3 x 16 ml)	HX083-YADEN (3 x 16 ml)	---	---	HX027-YCDEN (2 x 16 ml)	HX028-YCDEN (2 x 16 ml)	HX029-50DEN (4 x 13 ml)	HX010-07XEN (4ml)	HX615-YADEN (4no)	HX030-YADEN (3 x 16 ml)	HX035-YCDEN (15ml)	
QD410-YAXEN i6000 Diagnostics	200	---	---	HK111-20KN (20 ml)	HK083-20KN (20 ml)	---	---	HK518-20XN (20 ml)	HK519-20XN (20 ml)	HK520-10XN (5 x 10ml)	HK124-7KTN (1 x 4ml)	---	HK100-20KN (20ml)	---	
QD551-YCXEN NanoVIP	100	HX043-08XN (8.5 ml)	HX045-08XN (16ml)	HX047-08XN (3x8 ml)	HX049-08XN (3x 8 ml)	---	---	HX052-08XN (2x8 ml)	HX057-08XN (2x8 ml)	HX059-07XN (4 x 7 ml)	HX061-04XN (4ml)	HX616-10X (2no)	HX063-08XN (3 x 8ml)	---	
QD551-YCDEN NanoVIP	100	HX044-08XN (8.5 ml)	HX046-08XN (8.5 ml)	HX048-08XN (3 x 8 ml)	HX050-08XN (3 x 8 ml)	---	---	HX053-08XN (2 x 8 ml)	HX058-08XN (2 x 8 ml)	HX060-07XN (4 x 7 ml)	HX062-04XN (4ml)	HX617-10X (2no)	HX064-08XN (3 x 8 ml)	---	
QD552-YAXEN NanoVIP 300	100	HX031-12XN (12 ml)	HX032-12XN (12 ml)	HX026-20XN (1 x 20 ml)	HX083-20XN (1 x 20 ml)	---	---	HX027-20XN (1 x 20 ml)	HX028-20XN (1 x 20 ml)	HX029-20XN (2 x 20 ml)	HX010-02XN (2ml)	HY001-20X (3no)	HX030-20XN (1 x 20 ml)	---	

Note: It is recommended that reagents not be substituted across kit lot numbers.

Category	Detection Systems	Revision No.	F
Document No.	932-IHC2STN-Appendix	Release Date	17-Jun-2024

©2020, BioGenex Laboratories. All rights reserved.

 48810 Kato Road, Suite 100E & 200E, Fremont, CA 94538, Tel: +1 (800) 421-4149, Fax: +1 (510) 824-1490; support@biogenex.com, www.biogenex.com