

DATA SHEET

Xmatrix™ Coverslip

Doc. No.932-XT118, Rev. No. D
Date of Release: 11-Sep-2024

Intended Use

This product is intended for laboratory use only. BioGenex Coverslips are to be used on Xmatrix® Staining system.

Specifications

1. Coverslips are used on Xmatrix® Staining system after every time reagent is applied to optimize use and minimize wastage of micro-reagents.
2. Coverslips come in three configurations to accommodate to fit with different Micro Chamber slides provided by Biogenex: a single full-size test area of 25 x 40 mm, a single 2/3-size test area of 25 x 25 mm, and three 1/3-size test areas per slide, each compartment measuring 18 x 18 mm.

Items Supplied

One of the following:

Cat. Nos.	Description
XT118-50X	25 x 40 Coverslips (50 Nos)
XT118-YRK	25 x 40 Coverslips (500 Nos)

Applications

1. Coverslips are to be used with Xmatrix® Automated Staining system.
2. Coverslip application optimizes use and minimizes wastage of micro-reagents.
3. It also avoids the evaporation of micro-reagents during incubation period.