

DATA SHEET
17ml Empty Assay Vial (White) for Xmatrix

Doc. No.932-XT159-24X Rev. C
Date of Release: 11-Sep-2024

Intended Use

This product is intended for laboratory use only. The 17ml Empty Assay Vial (White) is for use on BioGenex Xmatrix Infinity Staining Systems.

Items Supplied:

Cat. No.	Description	Figure
XT159-24X	17ml Empty Assay Vial (White), 24/pk	
XT159-10X	17ml Empty Assay Vial (White), 10/pk	

Note:

1. Proper positioning of all vials in reagent racks is crucial since the pipette head is calibrated to this height and orientation.
2. Recalibration of the reagent vial position is required when using the vials for the first time.
3. When loading these vials into the Automated Staining Systems reagent racks, make sure that the vial neck is to the left.